


IKA Tube Mill 소용량 튜브 밀

Cat. No.	Model	Type	용량	EA Price ₩
M03-05-041	4180000	Tube Mill Control	40 ml	4,750,000
M03-05-031	20014733	Tube Mill 100 Control	100 ml	6,260,000

- Disposable Grinding Chamber
 - PP일회용 Chamber를 사용하여 사용이 편리함
 - Grinding 후 Sample 보관 및 이송용으로 사용 가능
 - No Clean Cost
- 본체에 덮개가 달려 있어 조용하고 안정함 • OLED Display, USB Interface
- IKA



Cat. No.	M03-05-041	M03-05-031
Type	Tube Mill Control	Tube Mill 100 Control
Motor Rating (in/output)	100 / 80 W	
Speed	5,000 ~ 25,000 rpm	
Useful Volume	40 ml	100 ml
Power On Time	3 min	
Overload Protection	Yes	
Circumferential Speed	65 m/s	
Max. Feed Hardness	5 Mohs	
Max. Feed Grain Size	10 mm	
Size / Weight	w180×d300×h170 mm / 2.7 kg	w180× d300× h212 mm / 2.8 kg
Power	230V, 50/60 Hz	

Disposable PP Grinding Transparent Chamber

Cat. No.	Model	Type	용량	PK Price ₩
M03-05-046	4425000	MT40.10	40 ml	10ea 210,000
M03-05-019	20003165	MMT40.1	40 ml, Reusable	1ea 440,000
M03-05-023	20008386	MT100.10	100 ml	10ea 240,000



IKA Analytical Mill 분석용 밀, A10 basic

Cat. No.	Model	Type	용량	Speed	EA Price ₩
M03-05-061	25001027	A10 basic	50 ml	25,000 rpm	4,190,000

- 구성: A10 basic, A10.1 Cutter, A10.4 Chamber Reduction
- For dry Grinding of Hard, Beittle, Soft and Fibrous Material
- Safety Device
 - 뚜껑이 닫혀 있어야 작동이 되며, 완전히 멈춘 후 뚜껑이 열림
 - Quick Stop 기능
- Intergrated Cooling Connection
- IKA

Cat. No.	M03-05-061
Type	A10 basic
Motor Rating (in/output)	300 / 240 W
Speed	25,000 rpm (fixed)
Useful Volume	50 ml
Duty Cycle On/Off	5 / 10 min
Overload Protection	Yes
Circumferential Speed	73 m/s
Max. Feed Hardness	5 Mohs
Max. Feed Grain Size	6 mm
Size / Weight	w130×d145×h250 mm / 3.13 kg
Power	220 ~ 240V, 50/60 Hz



Disposable PP Grinding Transparent Chamber

Cat. No.	Model	Type	Description	EA Price ₩
M03-05-582	25001680	A 10.1	Spare Cutter	120,000