

Vortex Mixer 보급형 볼텍스 믹서

| Cat. No. | Model | Speed | 작동 | EA Price ₩ |
|-------------|--------|---------------|-------------|------------|
| M04-238-521 | MX-S ① | 0 ~ 3,000 | 연속 or Touch | 185,100 |
| M04-238-528 | MX-E ② | 3,000 (Fixed) | Touch | 156,500 |

- 가장 보편적으로 사용되는 기본형 Vortex Mixer
- Orbital Motion
- Touch 또는 Continuous Operation
- 안정감 있는 작동을 위한 Silicon Feet 와 Steel Base
- VT1.1 Standard Attachment 포함 (for up to Ø30 mm Test Tube)
- DLAB

| Cat. No. | M04-238-521 | M04-238-528 |
|--------------------------|------------------------------|-------------------|
| Model | MX-S | MX-E |
| Motor Rating (In/Output) | 58w/10w | 12w/10w |
| Speed | 0 ~ 3,000 rpm | 3,000 rpm (Fixed) |
| Orbital Diameter | 4 mm / 4.8 mm | |
| Motor Type | Shaded-pole Motor | |
| Permissible On Time | 30 min | |
| Operation | Continuous / Touch Operation | Touch only |
| Size | w127×d160×h130 mm | w133×d133×h80 mm |
| Weight | 3.5 kg | 0.6 kg |
| Power | 220 V, 50/60 Hz | |

Optional Accessory

| Cat. No. | Model | Description | EA Price ₩ |
|-------------|-----------|--|------------|
| M04-238-911 | VT1.1 ③ | Standard Cup (기본 포함) | 8,600 |
| M04-238-915 | VT1.3 ④ | Universal Top Plate, Ø100 mm | 18,600 |
| M04-238-916 | VT1.3.6 | Pad for VT 1.3 Top Plate | 5,800 |
| M04-238-917 | VT1.3.1 ⑤ | Test Tube Adapter, Ø6 mm (48홀) | 12,800 |
| M04-238-919 | VT1.3.2 | Test Tube Adapter, Ø10 mm (15홀) | 12,800 |
| M04-238-923 | VT1.3.3 | Test Tube Adapter, Ø12 mm (12홀) | 12,800 |
| M04-238-927 | VT1.3.4 ⑥ | Test Tube Adapter, Ø16 mm (8홀) | 12,800 |
| M04-238-931 | VT1.3.5 | Test Tube Adapter, Ø20 mm (8홀) | 12,800 |
| M04-238-933 | VT1.2 | Tube Holding Rod, VT1 3.1 ~ VT1 ~ 3.5용 | 54,300 |

Microplate Mixer 보급형 마이크로 플레이트 믹서, MX-M

| Cat. No. | Model | Speed | 조절 | EA Price ₩ |
|-------------|-------|---------------|-------|------------|
| M04-238-531 | MX-M | 0 ~ 1,500 rpm | Scale | 355,800 |


- 컴팩트 타입의 96 Well Plate용 Shaker
- Continuous Operation, Orbital Motion
- PS1.2 Top Plate 기본 포함 • DLAB

| Cat. No. | M04-238-531 |
|----------------------------|-----------------------------------|
| Model | MX-M |
| Motor Rating (In / Output) | 18 / 10 W |
| Speed | 0 ~ 1,500 rpm |
| Orbital Diameter | 4.5 mm |
| Motor Type | External Rotor Brushless DC Motor |
| Speed Display | Scale |
| Operation | Continuous |
| Size | w260×d80×h150 mm |
| Weight | 3 kg |
| Power | 100 ~ 240V, 50/60 Hz |

Spare Accessory

| Cat. No. | Model | Description | EA Price ₩ |
|-------------|---------|-------------------------------------|------------|
| M04-238-966 | PS1.2 | Spare Top Plate, For 1 × Microplate | 8,400 |
| M04-238-963 | PS1.1 | Top Plate, Round, Ø100 mm | 8,400 |
| M04-238-969 | MXM-116 | Top Plate, For 2 × Microplate | 9,600 |

- PS1.1은 MX-S/F Vortex Mixer용 Test Tube Adapter와 같이 사용



- MXM-116 착용 시 0 ~ 1,000 rpm