

Vortex Mixer 볼텍스 믹서

Cat. No.	Model	Speed	EA Price ₩
M04-06-100	VM-10	Max. 3300 rpm	198,000

- Shaking 속도 조절 가능
- 고무 압착발을 적용하여 미끄러짐이나 진동을 최소화
- 3단계 구동 스위치. 정지, 연속동작 모드, Touch On 모드
- 직관적이고 편리한 Touch On 모드와 속도 조절 다이얼
- 2가지 Head Platform을 적용하여 다양한 형상과 용량의 Tube 사용 가능
- 소음이 적고 내구성이 좋은 Brushless Shade Motor 적용
- 구성 : Mixer, Rubber Platform Head, Pop-off Cup Head, Remove Applicator
- PL Insurance(제조물 배상 책임 보험)
- LK Labkorea

Cat. No.	M04-06-100		
Type	VM-10		
Motion	Orbital		
Speed Range	0 ~ 3,300 rpm		
Size	w131 x d185 x h125 mm (pop off cup head 제외)		
Material	Powder Coated Aluminum		
Operating Mode	Touch or Continuous		
Motor	Brushless Shade motor		
Max. Load	500 g		
Environmental Conditions	Temp. 5 ~ 40 °C, Humidity < 80 %		
Electric Supply	Power	1 Phase, 220Vac, 50 or 60Hz	
	Consumption	30 W	

Spare Accessory

Cat. No.	Model	Speed	EA Price ₩
M04-06-190	VM10RP	3 Inch Rubber Platform Head	12,000
M04-06-195	VM10PC	Pop-off Cup Head	10,000
M04-06-200	VM10RA	Remove Applicator	8,000



Vortex Mixer, IKA Lab dancer & Hula Dancer 볼텍스 믹서, Orbital Action

Cat. No.	Model	Type	Speed	EA Price ₩
M04-05-107	3365000	Lab dancer	2,800 rpm (Fixed)	340,000
M04-05-151	25004320	Hula Dancer Basic		740,000
M04-05-155	25003634	Hula Dancer Digital	1,200 ~2,800 rpm	870,000

- 작고 컴팩트한 사이즈
- 직경 30 mm 이내의 Test Tube, Centrifuge Tube, Eppendorf Vessel을 효율적으로 Shaking
- 탁월한 Orbital Mixing Action
- Base가 튼튼하게 제작되어 Shaking시 안전성이 탁월함
- IKA



Cat. No.	M04-05-107	M04-05-151	M04-05-155
Model	3365000	25004320	25003634
Type	Lab Dancer	Hula Dancer Basic	Hula Dancer Digital
Motor Rating (In/Output)	1.2 / 0.8 W	9.5 / 5 W	
Speed	2,800 rpm (fixed)		1,200 rpm~2,800 rpm
Shaking Stroke	4.5 mm		
Permissible Shaking Weight	0.1 Kg		
Protection Class	IP40	IP42	
Size/Weight	W100xD100xH70 mm / 0.7 kg	W123xD123xH42 mm / 0.7 kg	
Power	100 - 240 V, 50/60 Hz		