

## Vortex Mixer for Microtiter plate 볼텍스 믹서, IKA MTS 2/4

Cat. No.	Model	Type	Speed	EA Price ₩
M04-05-147	3208000	MTS 2/4	50 ~ 1100 rpm	3,080,000

- 2개 또는 4개의 Microtiter Plates 사용을 위한 특수한 Shaker
- Digital timer 기능으로 설정 시간의 동작이 끝나면 자동 알람
- Timer 기능은 1 ~ 99 분까지 설정이 가능
- Timer 기능을 사용하지 않을 경우 무한대로 사용가능함
- Electronic Speed Control
- Microtiter Plate용 Attachment 기본 포함
- Protection class acc. to DIN EN 60529 (IP21)
- IKA

Cat. No.	M04-05-147
Type	MTS 2/4
Motor Rating (In / Output)	35 / 13.5 W
Speed	50 ~ 1100 rpm
Orbital Diameter	3 mm
Speed display	Scale
Timer	Digital (Continuous / 1 ~ 99 min)
Operation	Continuous Operation & Timer
Size	w185×d320×h105 mm
Weight	3.2 kg
Power	220 ~ 240 V, 50/60 Hz



MTS 2/4

CE

## Vortex Mixer 볼텍스 믹서, IKA VXR basic

Cat. No.	Model	Type	Speed	EA Price ₩
M04-05-142	2819000	VXR basic	100 ~ 2200 rpm	2,280,000

- 넓은 범위의 Speed Range를 가진 Opto-electronically Controlled 되는 소형 Shaker
- 연속작동 모드에 적합
- 새로운 디자인과 개선된 Drive System
- 느린 Speed도 잘 조절됨
- Circular Shaking Motions (Orbital)
- Protection class acc. to DIN EN 60529 (IP21)
- IKA

Cat. No.	M04-05-136
Type	VXR basic
Motor Rating (In / Output)	35 / 13.5 W
Speed	100 ~ 2200 rpm
Orbital Diameter	4 mm
Speed display	Scale
Timer	None
Operation	Continuous Operation
Size	w157×d247×h130 mm
Weight	6.1 kg
Power	220 ~ 240 V, 50/60 Hz

### Optional Accessory

Cat. No.	Model	Type	Description	EA Price ₩
M04-05-632	0607200	VX 1 ①	One-hand Attachment	260,000
M04-05-636	0568900	VX 2 ②	Test Tube Attachment	1,130,000
M04-05-640	1618100	VX 2E ③	"Eppendorf" Attachment	900,000
M04-05-644	0953300	VX 7 ④	Dish Attachment	700,000
M04-05-648	0910400	VX 8 ⑤	Universal Attachment	1,670,000
M04-05-652	3627700	VX 11 ⑥	Basic Holder	370,000
M04-05-656	3659000	VX 11.1 ⑦	Test Tube Insert, Ø10 mm 70개	370,000
M04-05-660	3659100	VX 11.2	Test Tube Insert, Ø12 mm 41개	370,000
M04-05-664	3659200	VX 11.3	Test Tube Insert, Ø16 mm 32개	280,000
M04-05-668	3659300	VX 11.4 ⑧	Test Tube Insert, Ø20 mm 18개	280,000



VXR basic

CE

