

Plastic Mortar 플라스틱 유발

Cat. No.	Model	용량	Size	EA Price ₩
M05-59-201	110705	300 ml	od125×h80 mm	60,600
M05-59-203	110710	500 ml	od150×h90 mm	65,300

- 재질 : MF, 120°C 내열, Autoclavable
- Pestle(봉) 기본 포함
- Brand GmBH



Quartz Mortar 석영 유발

Cat. No.	Model	Top od x id	Height	EA Price ₩
M05-695-105	85-1511	od60 x id50 mm	20 mm	1,239,000
M05-695-110	85-1512	od70 x id55 mm	25 mm	1,544,600
M05-695-115	85-1513	od80 x id65 mm	30 mm	2,019,800
M05-695-120	85-1514	od90 x id70 mm	33 mm	2,630,800
M05-695-125	85-1515	od100 x id80 mm	35 mm	3,089,000
M05-695-130	85-1516	od120 x id100 mm	42 mm	5,567,100

- 재질 : 규소(SiO₂) 99.95% 이상
- Mortar 내부와 봉은 연이되어져 있음
- Sansyo



Stainless Steel Mortar And Pestle 스테인레스 유발 세트

Cat. No.	Model	용량	id	od	높이	EA Price ₩
M05-137-121	5400	165 ml	90 mm	100 mm	60 mm	155,500
M05-137-125	5401	330 ml	125 mm	135 mm	80 mm	189,500
M05-137-129	5402	450 ml	150 mm	160 mm	90 mm	225,200
M05-137-133	5403	1 L	180 mm	190 mm	105 mm	266,400

- 재질 : 18/10 Stainless Steel Antimagnetic, Chemical-Resistant Mortar,
18/9 Stainless Steel Pestle Head, Powder coated steel Handle
- Pestle(봉) 기본 포함
- Usbeck



Zirconia Mortar 지르코니아 유발

Cat. No.	Model	od	id	EA Price ₩
M05-196-232	37ZM-R705	60 mm	50 mm	1,098,000
M05-196-239	37ZM-R716	70 mm	60 mm	1,158,900
M05-196-246	37ZM-R727	80 mm	65 mm	1,326,600
M05-196-252	37ZM-R732	90 mm	70 mm	1,601,100
M05-196-259	37ZM-R734	100 mm	80 mm	1,906,000
M05-196-266	37ZM-R739	120 mm	105 mm	2,790,500
M05-196-273	37ZM-R748	150 mm	130 mm	3,476,600

- 94.8% ZrO₂ + 5.2% Y₂O₃
- Compressive Strength 5000 Mpa
- Material Density : 6.0 kg/m³
- Pestle(봉) 은 기본 포함
- LabCeramic®

