



## Silicon Oil - Silicon Fluid 실리콘 오일

Cat. No.	Model	내열 온도	점도(at 25°C)	포장	EA Price ₩
002-543-666	F200-10-1	-40 ~ 150°C	10 CS	1 kg	41,800
002-543-672	F200-10-4			4 kg	151,500
002-543-678	F200-10-16			16 kg	551,300
002-543-680	F200-50-4		50 CS	4 kg	131,300
002-543-682	F200-50-16			20 kg	476,100
002-543-684	F200-100-1		100 CS	1 kg	41,000
002-543-692	F200-100-4			4 kg	117,700
002-543-700	F200-100-16			20 kg	531,700

- Polydimethylsiloxane Polymer
- 무색 무향의 열매체유로 사용 가능한 실리콘 오일
- 휘발성이 낮아 오랫동안 사용 가능함 • Dow Corning



## High Temperature Silicon Oil 고온용 메틸 페닐 실리콘 오일

Cat. No.	Model	내열 온도	점도 (at 25°C)	포장	EA Price ₩
002-246-111	CD550-1	-40 ~ 232°C	125 CS	1 kg	116,800
002-246-118	CD550-4			4 kg	415,300

- “Phenylmethyl Siloxane”의 고온용 열매체
- 무색, 무취, 투명한 저점도 고온용 실리콘 오일
- Dow Corning



## Silicon Oil 실리콘 오일

- ShinEtsu

### 기본형 실리콘 오일

Cat. No.	Model	내열 온도	점도 (at 25°C)	포장	EA Price ₩
002-245-231	KF96-10	Max. 120°C	10 CS	1 Kg	55,100
002-245-240	KF96-50		50 CS		70,800
002-245-249	KF96-100	Max. 140°C	100 CS		21,900
002-245-258	KF96-1000		1000 CS		16,000

- 무색, 무취, 투명의 디메틸 폴리 실록산 구조를 가진 오일
- 열에 강한 실록산 결합과 유기질의 메틸기로 구성

### 내열용 실리콘 오일

Cat. No.	Model	내열 온도	점도 (at 25°C)	포장	EA Price ₩
002-245-311	KF54	-30 ~ 250°C	400 CS	1 kg	103,700

- 디메틸 실리콘 오일의 메틸기 일부를 페닐기로 치환한 메틸페닐계 실리콘 오일
- 오일베스등 주로 개방형 열매체유로 사용

### 내한용 실리콘 오일

Cat. No.	Model	내열 온도	점도 (at 25°C)	포장	EA Price ₩
002-245-411	KF50-100	-60 ~ 200°C	100 CS	1 kg	103,700

- 저온용 열매체유로 주로 사용됨



## Gold Vacuum Pump Oil, DirectTorr® 진공펌프 오일, Welch

Cat. No.	Model	용량	EA Price ₩
P06-750-995	8995G-11	1 L	140,000
P06-750-997	8995G-15	1 Gallon (3.8 L)	481,900

- Corrosive Gases or Vapors 용
- Welch Chenstar Pump와 같은 Belt-driven과 사용하기 알맞음
- Double Distilled Synthetic Hydrocarbon Oil
- No Additives or Inhibitors
- High Stability in Chemical Environments with Exceptionally Low Vapor Pressure
- Welch-ilmvac

