

pH / Ion / Conductivity Meter pH / Ion / Conductivity 미터

- 정밀형 pH, ORP, Ion, Conductivity Concentration Meter
- 빠르고 정밀하며 작동 및 거치가 편리함
- Large Color Display and Intuitive Touch Screen Operation
- ISM® (Intelligent Sensor Management) 기능
- 구성
 - 1) SD20/SD50 : Meter, Electrode Arm EasyPlace, Protect Cover, Manual, Quick Guide Power Supply, Declaration of Conformity and Test Certificate
 - 2) SD20/SD50-Kit : “SD20/SD50” 기본 구성품, InLab® Expert Pro-ISM 전극, pH 측정가이드, 1회용 Buffer 3개 (pH 4.01, 7.00, 9.21 and 10.00)
 - 3) 이외의 Kit : “SD20/SD50-Kit” 기본구성품, 사용 Kit에 따라 다른 전극으로 구성
 - 4) SD23 : Meter, Sensor Holder EasyPlace Triple, Protect Cover, Manual, Quick Guide Power Supply, Declaration of Conformity and Test Certificate
 - 5) SD23-Kit : “SD23” 기본 구성품, InLab® Expert Pro-ISM 전극, InLab® 731-ISM 전극, 1회용 Buffer 3개 (pH 4.01, 7.00, 9.21 and 10.00)
- Mettler Toledo / SevenCompact®



pH Meter

Cat. No.	Model	구분	전극	EA Price ₩
P08-162-201	SD20	Instrument	(포함 안됨)	1,916,800
P08-162-204	SD20-Kit	InLab® Expert Pro-ISM		2,278,600
P08-162-207	SD20-HA	InLab® Routine Pro-ISM		2,389,400
P08-162-210	SD20-Organic	InLab® Science Pro-ISM		2,727,600
P08-162-213	SD20-Pure H2O	InLab® Pure Pro-ISM		2,747,800
P08-162-216	SD20-Solids	InLab® Solids Pro-ISM		2,644,800
P08-162-219	SD20-Viscous	InLab® Viscous Pro-ISM		2,880,800
P08-162-222	SD20-LowVol	InLab® Micro Pro-ISM		2,799,000
P08-162-225	SD20-Stirrer	InLab® Expert Pro-ISM & EasyMix stir		2,778,000

pH / Ion Meter

Cat. No.	Model	구분	전극	EA Price ₩
P08-162-231	SD50	Instrument	(포함 안됨)	2,195,800
P08-162-234	SD50-Kit	InLab® Expert Pro-ISM		2,602,200
P08-162-237	SD50-HA	InLab® Routine Pro-ISM		2,714,000
P08-162-240	SD50-Organic	InLab® Science Pro-ISM		3,043,200
P08-162-243	SD50-Pure H2O	InLab® Pure Pro-ISM		3,208,800
P08-162-246	SD50-LowVol	InLab® Micro Pro-ISM		3,250,400
P08-162-249	SD50-F-Ion	PerfectION™ F 센서		4,032,000

pH / Conductivity Meter

Cat. No.	Model	구분	pH 전극	Conductivity 전극	EA Price ₩
P08-162-251	SD23	Meter Only	X	X	3,302,000
P08-162-254	SD23-Kit	Basic	InLab® Expert Pro-ISM	InLab® 731-ISM	4,495,400
P08-162-257	SD23-Pure H ₂ O	초순수용	InLab® Pure Pro-ISM	InLab® 741-ISM	4,870,400

• 자세한 사항은 280 Page 참조

SD20 / SD50	Range	Resolution	Accuracy
pH	-2.000 ~ 20.000	0.001 / 0.01 / 0.1	±0.002
mV	-2000.0 ~ 2000.0	0.1 / 1	±0.2
Ion Concentration(SD50)	0 ~ 999999 mg/L ppm	0.001 ~ 1	±0.5%
	0 ~ 100 mol/L	0.0001 ~ 1	±0.5%
	0 ~ 100.000 mmol/L	0.001 ~ 1	±0.5%
	0.000 ~ 100.000%	0.0001 ~ 1	±0.5%
	-2 ~ 20 pX	0.0001 ~ 1	±0.002
Conductivity Concentration(SD23)	0.001 μS/cm ~ 2000 mS/cm	0.001 ~ 1.000	±0.5%
Temp. MTC	-30 ~ 130°C	0.1°C	±0.1°C (0~100°C), 그외 ±0.3°C
Data Processing	Data Storage : 2000 Measurements Print Formats : Short, Normal, GLP		
Connector	pH/Ion (BNC), Temp. (Cinch/RCA), USB-A(프린터, 메모리 스틱), USB-B(PC), Mini-DIN(교반기)		
Size / Weight	205×195×65 mm / 850 g		

• Electrode Only, Cable은 별도 구매해야 함, Optional Accessory 참조

다음 페이지에서 계속...