

Portable Multi Meter 휴대용 멀티 미터, SevenGo Duo Pro™

- 최상의 요구를 만족하는 pH/ 이온/ 용존산소 미터 (IP67)
- 7개의 미리 설정된 보정점 중에서 5점 pH보정 또는 사용자가 정한 보정점에서 보정
- 500개의 GLP 데이터 저장 (날짜 / 시간, 센서ID, SN, 사용자 및 샘플 ID)
- 적외선 통신
- Intelligent Sensor Management (ISM®)
- Mettler Toledo

pH / 이온 / 용존 산소 미터, SG68

Cat. No.	Model	Standard Equipment	EA Price ₩
P02-162-250	SG68-B	Measuring Instrument	
P02-162-254	SG68-ELK	Electrode kit Inlab Expert Pro-ISM, Inlab 605-ISM(both IP67, 1.8m cable)	
P02-162-258	SG68-FK2	Electrode kit Inlab Expert Pro-ISM, Inlab 605-ISM(both IP67, 1.8m cable) and uGo™	
P02-162-262	SG68-FK5	Field kit with Inlab Expert Pro-ISM-5m, Inlab 605-ISM-5m(both IP67, 5m cable) and uGo™	

- pH, mV, 이온 농도, DO를 측정할 수 있는 휴대용 2 채널 미터

pH / 이온 / 전도도 미터, SG78

Cat. No.	Model	Standard Equipment	EA Price ₩
P02-162-282	SG78-B	Measuring Instrument	
P02-162-286	SG78-ELK	Electrode kit Inlab Expert Pro-ISM, Inlab 738-ISM(both IP67, 1.8m cable)	
P02-162-290	SG78-Fk2	Field kit with Inlab Expert Pro-ISM, Inlab 738-ISM(both IP 67, 1.8m cable) and uGo™	
P02-162-294	SG78-USP/EP	Field kit with Inlab Expert Pro-ISM, Inlab 738-ISM(both IP 67, 1.8m cable)	

- pH, mV, Conductivity, 전도도 TDS (Total Dissolved Solids), 염도 및 저항을 측정할 수 있는 휴대용 2 채널 미터

Model		SG68	SG78
pH	Range	-2.000 ~ 19.999	-2.000 ~ 19.999
	Resolution	0.001	0.001
	Accuracy	±0.002	±0.002
Temperature	Range	-5.0 ~ 130.0°C (ATC)	-30.0°C ~ 130.0°C (MTC)
	Resolution	0.1°C	0.1°C
	Accuracy	±0.2°C	±0.2°C
mV	Range	-1999 ~ 1999	-1999 ~ 1999
	Resolution	0.1	0.1
	Accuracy	±1	±0.1
Ion	Range	0.000 ~ 999.9%	-
	Resolution	3 digits	-
	Accuracy	±0.5%	-
Saturation	Range	0.0 ~ 600%	-
	Resolution	0.1 ~ 1	-
	Accuracy	±0.5%	-
Oxygen	Range	0.00 ~ 99.00	-
	Resolution	0.01	-
	Accuracy	±0.5%	-
Conductivity	Range	-	0.01 uS/cm ~ 1000 mS/cm
	Resolution	-	0.01
	Accuracy	-	±0.5%
TDS	Range	-	0.01 mg/L ~ 600 g/L
	Resolution	-	0.01 ~ 1
	Accuracy	-	±0.5%
Salinity	Range	-	0.00 ~ 80.00 ppt

공통사항

Sensor Inputs	BNC (10 ¹² Ω); NTC 30 kΩ (both IP67)
Out Puts	IR to Printer or PC via RS232 or USB (Only SG8)
Power Requirements	4×AA Batteries 1.5 V or NiMH Accumulators 1.3V
Operation Conitions	0°C ~ 40°C, 5%~80% rel. humidity (non-cond)
Package Size	w220×d90×h45 mm



Accessories for SG68 & SG78

Cat. No.	Model	Description	EA Price ₩
P02-162-411	51302332	Infrared-USB Adapter	
P02-162-415	51302333	Infrared-RS232 Adapter	

Other Accessory

Cat. No.	Model	Description	EA Price ₩
P02-162-455	51302325	Clip	88,600
P02-162-459	51302319	2-electrode Clip	94,800