

Ion Cell 이온 전극

• Mettler Toledo



Ion-Selective Electrode

Cat. No.	Model	Ion	측정범위	Membrane Type	pH Range	몸체 재질	Cable	EA Price ₩
P02-162-848	51107673	Lithium	10 ⁻⁵ ~ 10 ⁰ mol/L	Polymer	2 ~ 12 pH	POM/PVC	52300004	2,882,200
P02-162-852	51340900	Ammonium			2 ~ 9 pH			2,444,200
P02-162-856	51340500	Fluoride	10 ⁻⁶ ~ 10 ⁰ mol/L	Solid state	4 ~ 10 pH	POM		2,444,200
P02-162-860	51340263	Sodium	10 ⁻⁷ ~ 10 ⁰ mol/L	Na-glass	8 ~ 11 pH	Glass		1,172,400
P02-162-864	51107681	Cyanide	10 ⁻⁵ ~ 10 ⁰ mol/L	Solid state	4 ~ 13 pH	POM		2,882,200
P02-162-868	51107675	Sulfide	10 ⁻⁶ ~ 10 ⁰ mol/L					2,882,200
P02-162-872	51340400	Chloride	10 ⁻⁵ ~ 10 ⁰ mol/L	Polymer	2 ~ 13 pH	POM/PVC		2,356,800
P02-162-876	51340700	Potassium	10 ⁻⁶ ~ 10 ⁰ mol/L		2 ~ 12 pH			2,356,800
P02-162-880	51340600	Calcium		Premium 4-in-1	2 ~ 10 pH	POM		2,529,200
P02-162-884	51107870	Thiocyanate	10 ⁻⁵ ~ 10 ⁰ mol/L	Solid state	2 ~ 10 pH	POM		2,882,200
P02-162-888	51340800	Nitrate		Polymer	2 ~ 12 pH	POM/PVC	2,356,800	
P02-162-892	51107678	Copper	10 ⁻⁶ ~ 10 ⁰ mol/L	Solid state	2 ~ 8 pH	POM	2,882,200	
P02-162-896	51340300	Bromide			2 ~ 13 pH		2,529,200	
P02-162-900	51107676	Fluoroborate	10 ⁻⁵ ~ 10 ⁰ mol/L	Polymer	2 ~ 12 pH	Polymer	2,882,200	
P02-162-904	51107672	Cadmium			2 ~ 8 pH		2,882,200	
P02-162-908	51107680	Iodide	10 ⁻⁶ ~ 10 ⁰ mol/L	Solid state	1 ~ 13 pH	Solid state	2,882,200	
P02-162-912	51107674	Barium	10 ⁻⁵ ~ 10 ⁰ mol/L	Polymer	2 ~ 12 pH	Polymer	2,882,200	
P02-162-916	51107873	Lead			2 ~ 8 pH		2,882,200	

• 단독으로 사용할 수 없으며, 기준전극이 필요함. Spare Membrane Kit와 Filling Solution은 아래 참조

Spare Membrane Kit

Cat. No.	Model	Description	EA Price ₩
P02-162-935	51107687	Membrane Kit for Lithium	1,412,000
P02-162-938	51340012	Membrane Kit for Ammonium	1,569,800
P02-162-941	51340008	Membrane Kit for Fluoride	1,569,800
P02-162-944	51107695	Membrane Kit for Cyanide	1,412,000
P02-162-947	51107689	Membrane Kit for Sulfide	1,412,000
P02-162-950	51340007	Membrane Kit for Chloride	1,303,400
P02-162-953	51340010	Membrane Kit for Potassium	1,441,200
P02-162-956	51340009	Membrane Kit for Calcium	1,569,800
P02-162-959	51107871	Membrane Kit for Thiocyanate	1,399,600
P02-162-962	51340011	Membrane Kit for Nitrate	1,441,200
P02-162-965	51107692	Membrane Kit for Copper	1,412,000
P02-162-968	51340006	Membrane Kit for Bromide	1,569,800
P02-162-971	51107690	Membrane Kit for Fluoroborate	1,412,000
P02-162-974	51107686	Membrane Kit for Cadmium	1,412,000
P02-162-977	51107694	Membrane Kit for Iodide	1,412,000
P02-162-980	51107688	Membrane Kit for Barium	1,412,000
P02-162-983	51107874	Membrane Kit for Lead	1,399,600

Spare Filling Solution

Cat. No.	Model	Description	EA Price ₩
P02-162-989	51107881	Electrolyte for Lithium	122,200
P02-162-992	51340035	Electrolyte for Ammonium	74,400
P02-162-995	51340031	Electrolyte for Fluoride	74,400
P02-162-998	51107893	Electrolyte for Cyanide	122,200
P02-162-S04	51107894	Electrolyte for Sulfide	122,200
P02-162-S08	51340030	Electrolyte for Chloride	74,400
P02-162-S12	51340033	Electrolyte for Potassium	74,400
P02-162-S16	51340032	Electrolyte for Calcium	74,400
P02-162-S20	51107872	Electrolyte for Thiocyanate	121,200
P02-162-S24	51340034	Electrolyte for Nitrate	74,400
P02-162-S28	51107889	Electrolyte for Copper	122,200
P02-162-S32	51340029	Electrolyte for Bromide	74,400
P02-162-S36	51107890	Electrolyte for Fluoroborate	122,200
P02-162-S40	51107891	Electrolyte for Cadmium	122,200
P02-162-S44	51107898	Electrolyte for Iodide	122,200
P02-162-S48	51107892	Electrolyte for Barium	122,200
P02-162-S52	51107875	Electrolyte for Lead	121,200

pH Standard Solution pH 표준 용액

Cat. No.	Model	pH	Capacity	EA Price ₩
P02-162-666	51350002	2.00	250 ml Bottle	42,600
P02-162-669	51350004	4.01		40,400
P02-162-672	51350006	7.00		40,400
P02-162-675	51350008	9.21		40,400
P02-162-678	51350010	10.00		42,600
P02-162-681	51350012	11.00	Sachets 30×20 ml	48,400
P02-162-684	51302069	4.01		150,200
P02-162-687	51302047	7.00		150,200
P02-162-690	51302070	9.21		150,200

• 각 각의 버퍼는 COA (Certificate of Analysis) 와 같이 제공됨 (용기 라벨에 인쇄되어 있음)

• Mettler Toledo

