



## Disposable Dropping Pipet 일회용 드로핑 피펫 - 스포이드, Rayhope®

Cat. No.	Model	전체용량	Bulb	PK Price ₩	CS Price ₩
P03-35-066	DP-012	1.2 ml	0.9 ml	500ea 20,600	10pk 187,200
P03-35-072	DP-015*	1.5 ml	0.8 ml	500ea 20,600	10pk 187,200
P03-35-078	DP-030*	3.0 ml	2.1 ml	500ea 20,600	10pk 187,200
P03-35-084	DP-040	4.0 ml	3.1 ml	500ea 20,600	10pk 187,200
P03-35-090	DP-050*	5.0 ml	3.5 ml	500ea 20,600	10pk 187,200
P03-35-096	DP-055*	5.6 ml	1.9 ml	500ea 20,600	10pk 187,200
P03-35-102	DP-060	6.0 ml	3.0 ml	500ea 20,600	10pk 187,200
P03-35-108	DP-075*	7.5 ml	3.2 ml	500ea 20,600	10pk 187,200
P03-35-114	DP-230	23 ml	4.8 ml	100ea 20,600	10pk 187,200

- Inert, Nontoxic Polyethylene
- 유연성, 투명성, 내열성 (-50 ~ 90°C)이 뛰어남
- “\* ”는 눈금을 가지고 있음 (Graduated Pipet)
- Rayhope®

Model	전체길이 A	Bulb 길이 B	Bulb 외경 C	Stem 외경 D	눈금
DP-012	64 mm	27 mm	9.2 mm	2.5 mm	-
DP-015*	116 mm	26 mm	9.4 mm	3.2 mm	0.3 ml
DP-030*	140 mm	30 mm	13.2 mm	5.0 mm	1 ml
DP-040	155 mm	40.7 mm	13 mm	2.7 mm	-
DP-050*	150 mm	46.0 mm	13.2 mm	4.7 mm	1 ml
DP-055*	155 mm	35 mm	11 mm	6.3 mm	2 ml
DP-060	230 mm	45 mm	12 mm	4.5 mm	-
DP-075*	155 mm	42 mm	13 mm	7.5 mm	3 ml
DP-230	300 mm	67 mm	18 mm	8.5 mm	-



## Glass Dropping Spoid 유리 스포이드

Cat. No.	Model	용량	눈금	길이	PK Price ₩
P03-96-135	DS-G-03	3 ml	0.1 ml	160 mm	5ea 8,000
P03-96-142	DS-G-05	5 ml		180 mm	5ea 11,300
P03-96-149	DS-G-10	10 ml	0.2 ml	210 mm	5ea 13,300
P03-96-156	DS-G-20	20 ml		250 mm	5ea 25,100

- 눈금형 스포이드
- Rubber Bulb 포함
- LK Labkorea



## Disposable Micro Pipet 일회용 마이크로 유리 피펫

Cat. No.	Model	용량	Tol.	색상	PK Price ₩
P03-23-875	708707	1, 2, 3, 4, 5 ul	±0.3	White	250ea 55,500
P03-23-882	708709	10 ul		Orange	250ea 48,000
P03-23-889	708718	20 ul	±0.25	Black	250ea 48,000
P03-23-896	708722	25 ul		2 x White	250ea 48,000
P03-23-899	708728	20+40 ul		2 x Red	250ea 49,800
P03-23-903	9600144	44.7 ul		Violet	250ea 68,400
P03-23-910	708733	50 ul		Green	250ea 48,000
P03-23-917	708744	100 ul		Blue	250ea 49,800
P03-23-920	708745	50+100 ul			250ea 49,800
P03-23-924	708757	200 ul		Red	200ea 70,700
P03-23-931	9602144	44.7 ul		2 x Violet	250ea 80,700

- Complied with DIN ISO 7550
- Conformity Certified for Volume Greater than 5 ul
- ISO 규격에 따른 Color Code, with Ring Mark
- 길이 : 125 mm
- “9602144” 모델은 Heparine 처리된 제품임
- Brand GmBH

