

Standard Pasteur Pipet 표준 파스츄어 피펫, ISO7712

Cat. No.	Model	길이	Tip Length	PK Price ₩	CS Price ₩
P03-59-055	747715	150 mm	45 mm	250ea 20,700	4pk 77,600
P03-59-061	747720	230 mm	120 mm	250ea 22,600	4pk 84,600
P03-59-063	747725	270 mm	160 mm	250ea 54,800	4pk 204,600
P03-59-065	747730	310 mm	210 mm	250ea 64,600	4pk 241,200

- Sodalime Glass
- Capacity : Approximately 2 ml
- 노즐 내경 (Tip id) : $\varnothing 1.0$ mm, 상부 외경 : $\varnothing 7.0$ mm
- Suction Tube Length : 25mm
- Brand



Standard Pasteur Pipet 표준 파스츄어 피펫, ISO7712

- Sodalime Glass
- 노즐 외경 : $\varnothing 1.5 \pm 0.2$ mm
- 상부 외경 : $\varnothing 6.8 \sim 7.1$ mm, 두께 : $\varnothing 0.5 \sim 0.52$ mm
- Poulten & Graf

표준형

Cat. No.	Model	길이	PK Price ₩	CS Price ₩
P03-53-555	D810	150 mm	250ea 37,300	4pk 149,000
P03-53-561	D812	230 mm	250ea 41,000	4pk 164,000

면전 마개 부착형

Cat. No.	Model	길이	PK Price ₩	CS Price ₩
P03-53-588	D810/PP	150 mm	250ea 56,400	4pk 225,400
P03-53-594	D812/PP	230 mm	250ea 58,700	4pk 234,500



Disposable Borosilicate Pasteur Pipet, Kimble® 일회용 파스츄어 피펫, ASTM E732-80

Cat. No.	Model	길이	PK Price ₩	CS Price ₩
P03-137-121	63A53	146 mm (5.75")	250ea 33,700	4pk 125,800
P03-137-125	63B92	228.6 mm (9.00")	250ea 44,900	4pk 167,700

- 재질 : Borosilicate Glass 51 Expansion, USP Type I, ASTM E438
- Unplugged Type, 요청에 의해 Cotton-plugged Type 공급 가능
- DWK Lifescience / Kimble®



Polystyrene Pasteur Pipet, Sterile Polystyrene 멸균 파스츄어 피펫

Cat. No.	Model	길이	Packing	CS Price ₩
P03-132-012	3488-001	146 mm (5.75")	Individually Wrapped	200ea 197,200
P03-132-016	3488-003	228.6 mm (9")	Individually Wrapped	200ea 244,200
P03-132-020	3488-025	146 mm (5.75")	Bulk Packed	500ea 427,300
P03-132-024	3488-027	228.6 mm (9")	Bulk Packed	500ea 471,900

- Polystyrene 재질로 유리에 비해 안전성이 뛰어나
- Irradiation Sterilized, RNase-free, DNase-free and Non-pyrogenic
- Sterile : Clean Room Class 100,000의 환경에서 제조
- Non-Sterile, Cotton Plugged : Ready to Use
- Complies with USP Class VI
- Chemglass

