

Disposable Sterile Glass Pipet, Kimble® 일회용 멸균 유리 피펫

Cat. No.	Model	용량 (ml)	눈금 (ml)	Tol. (ml)	PK Price ₩	CS Price ₩
P03-137-209	72105-11100	1	0.1	±0.01	200ea 161,200 5pk	753,000
P03-137-213	72105-21100	2	0.01		125ea 115,500 4pk	431,600
P03-137-217	72105-5110	5	0.1	±0.1	100ea 134,400 5pk	627,800
P03-137-221	72105-10110	10			100ea 132,700 5pk	620,000

- 재질 : Borosilicate Glass 5.1, USP Type I, ASTM E438 Type I, Class B
- Sterile, Non-Pyrogenic, Plugged : Ready to Use
- Calibrated To Deliver (TD) to the Tip, Negative Graduation
- Designed from ASTM Specification E714
- Dispensing Shelf Pack 안에 개별 포장 되어 있음
- ASTM E380에 따른 Color Coding
- DWK Lifescience / Kimble®



Disposable Serological Glass Pipet, Kimble® 일회용 유리 피펫, 비멸균

Cat. No.	Model	용량 (ml)	눈금 (ml)	Tol. (ml)	PK Price ₩	CS Price ₩
P03-137-173	72120-11100	1 ml	0.01	±0.03	250ea 186,200 4pk	695,800
P03-137-177	72120-21100	2 ml		±0.06	250ea 199,500 2pk	372,800
P03-137-181	72120-5110	5 ml	0.1	±0.15	100ea 109,200 5pk	510,100
P03-137-185	72120-10110	10 ml		±0.30	100ea 117,600 5pk	549,400

- 재질 : Borosilicate Glass 5.1, USP Type I, ASTM E438 Type I, Class B
- Non-Plugged, Non-Sterile
- Calibrated to Deliver (TD) to the Tip, Negative Graduation
- 혈청용 피펫은 마지막 방울을 붙여서 날려버려야 함
- ASTM E1380에 따른 Color Coding
- DWK Lifescience / Kimble®



Disposable PS Serological Pipet 일회용 플라스틱 피펫

Cat. No.	Model	용량	눈금	Color Code	PK Price ₩
P03-61-001	LK-H7201Y	1ml	0.01 ml	Yellow	1000ea 195,900
P03-61-003	LK-H7202G	2ml		Green	800ea 194,400
P03-61-005	LK-H7205B	5ml	0.1 ml	Blue	300ea 85,100
P03-61-010	LK-H7210R	10ml		Red	200ea 62,300
P03-61-025	LK-H7225P	25ml	0.2 ml	Purple	200ea 101,000
P03-61-050	LK-H7250B	50ml	1 ml	Black	100ea 109,000

- 높은 투명도를 제공하기 위해 polystyrene로 제작
- 6가지 색상의 Color Code, Individual Wrapper (개별 멸균포장)
- 흡입부 내부에 필터 삽입된 제품
- Non-cytotoxic, Non-pyrogenic
- 선명한 Ascending & Descending Scales 인쇄
- LK Labkorea



Disposable Serological Pipet 일회용 플라스틱 피펫, Sterile PS

Cat. No.	Model	용량	눈금	Color Code	CS Price ₩
P03-608-010	83.1251.001	1ml	0.01 ml	Yellow	1000ea 150,000
P03-608-013	83.1252.001	2ml		Green	1000ea 150,000
P03-608-016	83.1253.001	5ml	0.1 ml	Blue	500ea 110,000
P03-608-019	83.1254.001	10ml		Orange	500ea 110,000
P03-608-022	83.1685.001	25ml	0.2 ml	Red	25ea 80,000
P03-608-025	83.1256.001	50ml	0.5 ml	Purple	30ea 80,000

- 재질 : PS
- 멸균 개별 포장, Non-cytotoxic, Non-pyrogenic
- Antistatic Package
- Sarstedt

