

Graduated Pipet, Total Delivery, Type 2 메스 전량 피펫, Class AS + Batch 보증서

- AR-GLAS[®], Type 2, Serological Pipet
- Complied with DIN EN ISO 835 (Waiting Time : 5초)
- 하부 Tip까지의 용량이 상부에 표시

Class AS, with Batch Certificate / 보증서 포함, BLAUBRAND[®]

Cat. No.	Model	용량	눈금	Tol. (ml)	Color Code ISO1769	길이 (mm)	Cotten Plug	EA Price ₩	PK Price ₩
P03-59-601	27816	0.5 ml	0.01 ml	0.006	2×Yellow	360	No	7,800	12ea 87,300
P03-59-605	27817	1 ml		0.007	Yellow	360	No	7,100	12ea 79,300
P03-59-609	27819	2 ml	0.02 ml	0.01	Black	360	No	7,100	12ea 79,300
P03-59-613	27822	5 ml	0.1 ml	0.03	Blue	360	Yes	8,100	12ea 90,200
P03-59-617	27823	10 ml		0.05	Orange	360	Yes	8,400	12ea 93,600
P03-59-621	27824	20 ml		0.1	2×Yellow	360	Yes	15,700	6ea 87,800
P03-59-625	27825	25 ml	White		450	Yes	15,700	6ea 87,800	
P03-59-629	27826	50 ml	0.5 ml	0.2	-	450	Yes	28,700	6ea 160,700

- Conformity Certified, Calibrated to Deliver (TD.Ex)
- 보증서(Batch Certificate) 포함, Batch Number가 Pipet에 Mark되어짐
- Brand GmbH

Class AS, without Batch Certificate / 보증서 불포함

Cat. No.	Model	용량	눈금	색상 코드	EA Price ₩	PK Price ₩
P03-74-321	434116314	1 ml	0.01 ml	Yellow	3,000	10ea 27,100
P03-74-325	434116516	2 ml	0.02 ml	Black	3,000	10ea 27,100
P03-74-329	434116618	5 ml	0.05 ml	Red	3,000	10ea 27,100
P03-74-333	434116719	10 ml	0.1 ml	Orange	3,400	10ea 30,100
P03-74-337	632434116722	20 ml		2×Yellow	6,000	5ea 27,100
P03-74-341	434116723	25 ml		White	6,000	5ea 27,100
P03-74-345	434116825	50 ml	0.2 ml	Black	14,800	5ea 67,200
P03-74-349	434116831	100 ml	0.2 ml	Red	29,100	5ea 132,000

- 보증서가 포함되지 않았으나, Class AS의 정밀도를 가지는 보급형 제품임
- Techno Glass

Graduated Pipet, Total Delivery, Type 2, USP Grade, BLAUBRAND[®] 메스 전량 피펫, Class AS + USP Batch 보증서

Cat. No.	Model	용량 (ml)	눈금 (ml)	Tol. (ml)	Color Code ISO1769	길이 (mm)	Cotton Plug	EA Price ₩	PK Price ₩
P03-59-830	27506	1	0.01	0.007	Yellow	360	Yes	8,700	12ea 97,100
P03-59-834	27509	2	0.02	0.01	Green		Yes	8,700	12ea 97,100
P03-59-838	27511	5	0.05	0.02	Red		Yes	10,200	12ea 114,200
P03-59-842	27513	10	0.1	0.03	Orange		Yes	10,200	12ea 114,200

- AR-GLAS[®], Type 2, Serological Pipet
- Complied with DIN EN ISO 835 (Waiting Time : 5초)
- 하부 Tip까지의 용량이 상부에 표시
- Conformity Certified, Calibrated to Deliver (TD.Ex)
- USP Batch Certificate 포함
- Brand GmbH

ASTM Serialized & Certified TD Serological Pipet, Kimble[®] ASTM 메스 전량 피펫, Class A + 개별 보증서

Cat. No.	Model	용량	눈금	Tol.	Color Code	PK Price ₩	CS Price ₩	
P03-137-225	37025-11100	1 ml	0.01 ml	0.01 ml	Yellow	6ea 366,800	2pk 685,500	
P03-137-229	37025-2	2 ml			Green	6ea 372,700	2pk 696,600	
P03-137-233	37025-5	5 ml	0.1 ml	0.02 ml	Blue	6ea 330,200	2pk 617,200	
P03-137-237	37025-10	10 ml			0.03 ml	Orange	6ea 442,600	2pk 827,200
P03-137-241	37025-25	25 ml			0.05 ml	White	1ea 129,600	6pk 726,300

- KIMAX[®] Borosilicate Glass 3.3, USP Type I, Class A
- ASTM Specification E1293, Style I, Class A
- Brown Graduation
- 개별 피펫에 Serial No.가 Mark되어 지며, USP 개별 보증서 포함
- DWK Lifescience / Kimble[®]

