

Volumetric Pipet, Bulb Pipet, One Mark 홀 피펫, Class AS

- AR-GLAS[®], 1 Mark, Class AS
- Complied with DIN EN ISO 648 (Waiting Time : 5초)

Class AS, with Batch Certificate / 보증서 포함, BLAUBRAND[®]

Cat. No.	Model	용량	Tol.	Color Code	길이	EA Price ₩	PK Price ₩
P03-59-501	29701	0.5 ml	±0.005 ml	2 X Black	300 mm	5,700	12ea 63,100
P03-59-503	29702	1 ml	±0.008 ml	Blue		6,000	12ea 66,600
P03-59-505	29703	2 ml	±0.010 ml	Orange	330 mm	6,500	12ea 72,500
P03-59-507	29704	2.5 ml		-		7,400	6ea 41,300
P03-59-509	29705	3 ml	±0.015 ml	Black	400 mm	6,400	6ea 35,500
P03-59-511	29706	4 ml		2 X Red		6,500	6ea 36,100
P03-59-513	29707	5 ml		White		6,600	6ea 37,000
P03-59-515	29708	6 ml		2 X Orange		9,500	6ea 52,800
P03-59-517	29709	7 ml	±0.02 ml	2 X Green	440 mm	9,500	6ea 52,800
P03-59-519	29710	8 ml		Blue		9,500	6ea 52,800
P03-59-521	29711	9 ml	±0.03 ml	Black	510 mm	9,500	6ea 52,800
P03-59-523	29712	10 ml		Red		8,300	6ea 46,400
P03-59-525	29713	15 ml		Green		10,700	6ea 59,900
P03-59-527	29714	20 ml	±0.05 ml	Yellow	520 mm	10,200	6ea 57,000
P03-59-529	29715	25 ml		Blue		11,200	6ea 62,500
P03-59-531	29716	30 ml	±0.08 ml	Black	540 mm	11,900	6ea 66,500
P03-59-533	29717	40 ml		White		13,800	6ea 77,100
P03-59-535	29718	50 ml	±0.10 ml	Red	585 mm	12,600	6ea 70,700
P03-59-537	29719	100 ml		Yellow		18,100	6ea 101,200

- Conformity Certified, Calibrated to Deliver(TD, Ex)
- 보증서(Batch Certificate) 포함, Batch Number가 Pipet에 Mark되어짐
- Brand GmbH

Class AS, without Batch Certificate / 보증서 불포함

Cat. No.	Model	용량	눈금	Color Code	EA Price ₩	PK Price ₩
P03-74-533	433111014	1 ml	0.008 ml	Blue	4,600	10ea 41,000
P03-74-537	433111016	2 ml	0.010 ml	Orange	4,600	10ea 41,000
P03-74-548	433111018	5 ml	0.015 ml	White	4,800	10ea 43,600
P03-74-551	433111019	10 ml	0.020ml	Red	5,900	10ea 53,600
P03-74-554	433111021	15 ml	0.030 ml	Green	7,300	5ea 33,000
P03-74-557	433111022	20 ml		Yellow	8,100	5ea 36,600
P03-74-561	433111023	25 ml		Blue	8,200	5ea 36,900
P03-74-569	433111025	50 ml	0.050 ml	Red	11,600	5ea 52,500
P03-74-573	433111030	100 ml	0.080 ml	Yellow	21,900	5ea 99,500

- 보증서가 포함되지 않았으나, Class AS의 정밀도를 가지는 보급형 제품임
- Techno Glass

PP Volumetric Pipet 플라스틱 홀 피펫, Class B

Cat. No.	Model	용량	정밀도	길이	EA Price ₩	PK Price ₩
P03-20-601	30002	1 ml	±0.02 ml	300 mm	19,300	12ea 224,100
P03-20-606	30003	2 ml			19,700	12ea 229,600
P03-20-611	30007	5 ml	±0.03 ml	440 mm	36,100	6ea 209,900
P03-20-616	30012	10 ml	±0.04 ml		38,000	6ea 221,300
P03-20-621	30015	25 ml	±0.06 ml	450 mm	39,800	6ea 231,400
P03-20-626	30018	50 ml	±0.10 ml	460 mm	41,000	6ea 238,400

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- Calibrated to meet class B Tolerances
- 60°C 이상의 온도에 노출되면 정확성에 영향을 줄 수도 있음
- Brand GmbH

