

Volumetric Pipet, Bulb Pipet, One Mark

홀 피펫, Class AS + USP Batch 보증서

Cat. No.	Model	용량	Tol.	Color Code	길이	EA Price ₩	PK Price ₩
P03-59-539	30601	0.5 ml	±0.005 ml	2 X Black	300 mm	13,400	12ea 149,300
P03-59-541	30602	1 ml	±0.006 ml	Blue	330 mm	10,500	12ea 117,700
P03-59-543	30603	2 ml		Orange		10,500	12ea 117,700
P03-59-545	30605	3 ml	±0.01 ml	Black	400 mm	12,900	6ea 72,200
P03-59-547	30606	4 ml		2 X Red		14,900	6ea 83,500
P03-59-549	30607	5 ml		White		12,900	6ea 72,200
P03-59-551	30608	6 ml		2 X Orange		17,600	6ea 98,600
P03-59-553	30609	7 ml	±0.02 ml	2 X Green	440 mm	17,600	6ea 98,600
P03-59-555	30610	8 ml		Blue		17,600	6ea 98,600
P03-59-557	30611	9 ml		Black		17,600	6ea 98,600
P03-59-559	30612	10 ml	±0.03 ml	Red	510 mm	13,300	6ea 74,100
P03-59-561	30613	15 ml		Green		16,400	6ea 91,700
P03-59-563	30614	20 ml		Yellow		16,600	6ea 92,700
P03-59-565	30615	25 ml		Blue		16,900	6ea 94,800
P03-59-567	30616	30 ml	±0.05 ml	Black	520 mm	19,700	6ea 110,100
P03-59-569	30617	40 ml		White		24,200	6ea 135,500
P03-59-571	30618	50 ml		Red		22,300	6ea 124,600
P03-59-574	30619	100 ml	±0.08 ml	Yellow	585 mm	30,700	6ea 172,000

- AR-GLAS®, 1 Mark, Class AS
- Complied with DIN EN ISO 648 (Waiting Time : 5초)
- USP Batch Certificate
- Brand GmbH



ASTM Serialized & Certified Volumetric Pipet

ASTM 홀 피펫, Class A + 개별 보증서, Kimble®

Cat. No.	Model	용량	Tol.	길이	Color Code	PK Price ₩
P03-137-245	37006-510	0.5 ml	0.006 ml	310 mm	2 x Black	6ea 187,500
P03-137-249	37006-1	1 ml			Blue	6ea 187,300
P03-137-253	37006-1510	1.5 ml		335 mm	-	6ea 226,600
P03-137-257	37006-2	2 ml			Orange	6ea 240,500
P03-137-261	37006-2510	2.5 ml	0.01 ml	350 mm	-	6ea 256,900
P03-137-265	37006-3	3 ml			Black	6ea 278,400
P03-137-269	37006-4	4 ml		2 x Red	6ea 305,100	
P03-137-273	37006-5	5 ml		365 mm	White	6ea 314,300
P03-137-277	37006-6	6 ml			2 x Orange	6ea 315,500
P03-137-281	37006-7	7 ml		400 mm	2 x Green	6ea 320,600
P03-137-285	37006-8	8 ml	Blue		6ea 325,100	
P03-137-289	37006-9	9 ml	Black		6ea 326,500	
P03-137-293	37006-10	10 ml	0.02 ml	460 mm	Red	6ea 327,300
P03-137-297	37006-12	12 ml			-	3ea 174,000
P03-137-301	37006-15	15 ml	0.03 ml	455 mm	Green	3ea 184,300
P03-137-305	37006-20	20 ml			-	3ea 207,800
P03-137-309	37006-25	25 ml		505 mm	Blue	3ea 207,800
P03-137-313	37006-30	30 ml			Black	3ea 215,200
P03-137-317	37006-40	40 ml	0.05 ml	525 mm	White	3ea 242,900
P03-137-321	37006-50	50 ml			Red	3ea 247,000
P03-137-325	37006-75	75 ml	0.08 ml	550 mm	Green	3ea 285,500
P03-137-329	37006-100	100 ml			Yellow	3ea 298,800

- KIMAX® Borosilicate Glass 3.3, USP Type I, Class A
- ASTM Specification E969, Class A (75 & 200 ml : ASTM E542)
- Brown Graduation
- 개별 피펫에 Serial No.가 Mark되어 지며, USP Individual Certificate 포함
- DWK Lifescience / Kimble®

