

Pipet Bulb, Spoid Bulb 피펫 벌브


- 스포이드 및 드로핑 바틀용으로 적합
- 실리콘, 라텍스 재질
- LK Labkorea

Silicon Spoid Bulb

Cat. No.	Model	용량	Size (Ø × h)	PK Price ₩	CS Price ₩
P03-11-320	SB-S02	2 ml	16×41 mm	10ea 4,800	5pk 21,800
P03-11-326	SB-S03	3 ml	24×53 mm	10ea 5,500	5pk 26,100
P03-11-332	SB-S05	5 ml	26×56 mm	10ea 6,100	5pk 29,000
P03-11-338	SB-S10	10ml	28×58 mm	10ea 8,500	5pk 40,400

Large Spoid Bulb, Rubber

Cat. No.	Model	용량	Size (Ø × h)	EA Price ₩
P03-11-360	SB-R25	25 ml	43×62 mm	1,500
P03-11-365	SB-R75	75 ml	64×81 mm	2,500



PVC Spoid PVC 스포이드, Pasteur Pipet용

Cat. No.	Model	Size	PK Price ₩
P03-53-955	D813	Ø12×L44 mm	100ea 66,700


- Ribbed PVC
- LK Labkorea



Pipet Filler 피펫 필러

Cat. No.	Model	Type	용량	Valve 수	EA Price ₩
P03-59-701	25300	Standard	up to 10 ml	3	17,300
P03-59-703	25315	Universal	up to 100 ml	3	17,800
P03-59-705	25400	Flip	up to 100 ml	2	21,400


- 재질 : Natural Rubber
- 3 Valves Type : Filling (S), Delivery (E), Release Air Valve (A)
- 2 Valves Type : Filling (↑), Delivery (↓)
- Release Air Valve 없이 공기 배출이 가능한 Type
- Brand GmbH



Pipet Filler 피펫 필러

Cat. No.	Model	Type	EA Price ₩
P03-11-452	PF-D01	국산 ①	3,900
P03-20-911	PWY145	외산 Type I ②	7,500

- 거의 모든 규격의 피펫에 공용으로 사용되는 피펫 보조기구
- 흡입/분주 조절식으로 Natural Rubber로 만들어짐
- LK Labkorea



Pipet Pump, Pi-Pump 피펫 펌프

Cat. No.	Model	용량	색상	EA Price ₩
P03-11-555	108002	2 ml	Blue	12,700
P03-11-561	108010	10 ml	Green	12,700
P03-11-567	108025	25 ml	Red	14,800

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- Delivery Valve Lever와 Thumbwheel로 조절
- 흡입/분주 조절이 자유롭고 용이함
- Rayhope®

