

"Brand" Multi-Channel Micropipet 멀티 채널 마이크로 피펫, Transferpette® multi channel

- Microplate Format을 위한 최상의 피펫임
- 구성 : Transferpette® S -8 or -12, Conformity Certified, Performance Certificate, 2 Tip-Box filled with Pipette Tips, 1 Shelf/rack Mount, 1 Reagent Reservoir, Mounting Tool, Labeling Foil, Silicone Grease, 1 Set of Sealing Ring
- CE-IVD-Compliant
- Tip 선택시 "Brand" Tip Selection Guide (773 Page) 참조
- 기타 자세한 사항은 Transferpette® S와 동일함
- Brand GmBH



Transferpette® S-8

Cat. No.	Model	Type	용량	Accuracy	EA Price ₩
P09-59-223	705900	M8-10	0.5 ~ 10 ul	± 1.6%	1,328,300
P09-59-227	705906	M8-50	5 ~ 50 ul	± 0.8%	1,328,300
P09-59-231	705908	M8-100	10 ~ 100 ul		1,328,300
P09-59-235	705910	M8-200	20 ~ 200 ul		1,328,300
P09-59-239	705912	M8-300	30 ~ 300 ul	± 0.6%	1,328,300

Transferpette® S-12

Cat. No.	Model	Type	용량	Accuracy	EA Price ₩
P09-59-263	705920	M12-10	0.5 ~ 10 ul	± 1.6%	1,522,200
P09-59-267	705926	M12-50	5 ~ 50 ul	± 0.8%	1,522,200
P09-59-271	705928	M12-100	10 ~ 100 ul		1,522,200
P09-59-275	705930	M12-200	20 ~ 200 ul		1,522,200
P09-59-279	705932	M12-300	30 ~ 300 ul	± 0.6%	1,522,200

Rack for Transferpette® S & Transferpette® S-8 & -12

Cat. No.	Model	Description	EA Price ₩
P09-59-292	704807	Bench Top Rack ①	136,500
P09-59-907	704811	Shelf / Rack Mount ②	22,000
P09-59-909	704812	Shelf / Wall Mount ③	39,200
P09-59-297	703440	Individual Stand for Multi-Channel ④	59,500



- Mechanical Single & Multi-Channel Pipetter용 : ①, ②, ③
- Multi-Channel Pipetter용 : ④

"Brand" Tip Selection Guide

Pipette tips Filter tips Volume range	Transferpette® single channel Nominal volume													Transferpette® multichannel Nominal volume												
	1 µl	2.5 µl	5 µl	10 µl	10 µl	20 µl	20 µl	25 µl	50 µl	100 µl	200 µl	250 µl	500 µl	1000 µl	2 ml	5 ml	10 ml	10 µl	20 µl	25 µl	50 µl	100 µl	200 µl	300 µl		
0.1 - 20 µl	✓	✓ ₂₎	✓ ₂₎															✓	✓							
0.5 - 20 µl		✓	✓			✓												✓	✓							
1 - 50 µl		✓ ₂₎	✓			✓												✓	✓							
2 - 200 µl ¹⁾			✓			✓		✓	✓	✓	✓									✓	✓	✓	✓	✓	✓	✓
5 - 300 µl			✓			✓		✓	✓	✓	✓									✓	✓	✓	✓	✓	✓	✓
50 - 1000 µl ¹⁾												✓	✓	✓												
0.5 - 5 ml												✓	✓	✓												
1 - 10 ml																	✓									
0.1 - 1 µl	✓	✓ ₂₎	✓ ₂₎															✓	✓							
0.5 - 10 µl		✓	✓			✓												✓	✓							
1 - 20 µl			✓			✓												✓	✓							
2 - 20 µl			✓			✓		✓	✓	✓	✓									✓	✓	✓	✓	✓	✓	✓
5 - 100 µl			✓			✓		✓	✓	✓	✓									✓	✓	✓	✓	✓	✓	✓
5 - 200 µl											✓									✓	✓	✓	✓	✓	✓	✓
50 - 1000 µl													✓	✓												

✓ = Tip volume less than pipette's nominal volume
 1) Bulk tips are yellow or blue colored, racked tips are colorless in a yellow-colored or blue-colored tip tray
 2) Calibration and possible adjustment needed