

"Brand" Electronic Pipette 전동식 마이크로 피펫

Cat. No.	Model	용량	Subdivision	Accuracy	EA Price ₩
P09-59-744	705299	0.5 ~ 10 ul	0.01 ul	±1.0%	768,000
P09-59-748	705300	2 ~ 20 ul	0.02 ul	±1.0%	768,000
P09-59-752	705303	20 ~ 200 ul	0.2 ul	±0.8%	768,000
P09-59-756	705306	50 ~ 1000 ul	1.0 ul	±0.6%	768,000
P09-59-760	705307	500 ~ 5000 ul	5.0 ul	±0.6%	768,000

- 인체공학적 설계로 편리하고 쉬운 조작 (TUV 인증)
- Transferpette S와 같은 전량으로 사용이 편리함
- 간단하고 쉬운 Menu로 조작이 간편함
- Finger Rest가 조절 가능하여 손의 피로를 줄일수 있음
- Corrosion-resistant Piston and Ejector
- Five Convenient Program
 - Pipetting (PIP Mode)
 - Mixing of Sample (PIPMix Mode) : 반복적인 Pipetting, Dispensing
 - Reverse Pipetting (revPIP Mode) : High Viscosity or High Vapor Pressure Liquid용
 - Pipetting with Electrophoresis (GEL Mode) : For Loading of Gels
 - Dispensing (DISP Mode)
- Easy Calibration으로 자체 Calibration
- Battery Refresh Function, 4000 cycle pipetting (NiMH Battery)
- Pipette Shafts / Manifold re Completely Autoclavable at 121°C (20 min)
- Vapor pressure : up to 500 mbar, Viscosity: 260 mPa s
- 구성 : Transferpette® Electronic, Conformity Certified, Performance Certificate, Battery, Power Supply Unit, Silicon Oil
- Power : 110~240V, 50~60Hz
- Brand GmBH

Accessory

Cat. No.	Model	Description	EA Price ₩
P09-59-781	705390	3-Device Charging Stand ①	464,300
P09-59-784	705385	Individual Stand, up to 1000 ul	117,600
P09-59-787	705386	Individual Stand, 500 ~ 5000 ul ②	117,600
P09-59-790	705382	Shelf / Wall Mount	58,200
P09-59-793	705383	Shelf / Rack Mount	36,300



"Brand" Electronic Multi Channel Pipette 전동식 멀티 채널 마이크로 피펫

- 기타 자세한 사항은 Transferpette® Electronic의 특징 및 기능과 동일함
- Manifold Rotated by 360°
- 구성 : Transferpette® -8/-12 Electronic, Battery, Power Supply Unit, Device Stand, Tip-Box SL, Refill Unit, Conformity Certified, Performance Certificate, Reagent Reservoir, Silicon Oil
- Brand GmBH

Transferpette®-8 Electronic ①

Cat. No.	Model	용량	Subdivision	Accuracy	EA Price ₩
P09-59-822	705399	0.5 ~ 10 ul	0.01 ul	±1.2%	1,551,200
P09-59-826	705400	1 ~ 20 ul	0.02 ul	±1.0%	1,551,200
P09-59-830	705403	5 ~ 100 ul	0.1 ul	±0.8%	1,551,200
P09-59-834	705404	10 ~ 200 ul	0.2 ul	±0.8%	1,551,200
P09-59-838	705406	15 ~ 300 ul	0.5 ul	±0.6%	1,551,200

Transferpette®-12 Electronic ②

Cat. No.	Model	용량	Subdivision	Accuracy	EA Price ₩
P09-59-862	705449	0.5 ~ 10 ul	0.01 ul	±1.2%	1,792,300
P09-59-866	705450	1 ~ 20 ul	0.02 ul	±1.0%	1,792,300
P09-59-870	705453	5 ~ 100 ul	0.1 ul	±0.8%	1,792,300
P09-59-874	705454	10 ~ 200 ul	0.2 ul	±0.8%	1,792,300
P09-59-878	705456	15 ~ 300 ul	0.5 ul	±0.6%	1,792,300



Multi Channel Pipette Stand "703440"은 773 Page 참조

