

## "Brand" Piston-Operated Pipette 피스톤식 피펫, 고정도

Cat. No.	Model	용량	눈금	정밀도	EA Price ₩
P09-59-101	701807	2.5 ~ 10 ul	0.01 ul	±1.0%	746,500
P09-59-103	701812	5 ~ 25 ul	0.1 ul	±0.8%	746,500
P09-59-105	701817	10 ~ 50 ul			746,500
P09-59-107	701822	20 ~ 100 ul			746,500
P09-59-109	702804	100 ~ 500 ul	1.0 ul	±0.5%	746,500
P09-59-111	702806	200 ~ 1000 ul			746,500
P09-59-113	702810	1000 ~ 5000 ul			746,500
P09-59-115	702812	2000 ~ 10000 ul	10.0 ul		746,500

- 용도 :
  - Media tend to Foam : Surface Solution
  - Media with High Vapor Pressure : Alcohol, Ether, Hydrocarbons
  - Highly Viscous Media & Media with High Density
- 사용조건 :
  - Density : Up To 13.6 g/cm<sup>3</sup>
  - Viscosity : Up To 140,000 mm<sup>2</sup>/s
  - Vapor Pressure : Up To 500 mbar
- Working Temperature Range : 15°C ~ 40°C
- Conformity Certified
- 구성 : Transfepettor, Conformity Certified, Performance Certificate
- Brand GmBH

### Accessory

#### Cap

Cat. No.	Model	사용 용량	재질	PK Price ₩
P09-59-121	701902	for 10 ul	Glass ①	100ea 74,600
P09-59-123	701906	for 25 ul		100ea 74,600
P09-59-125	701908	for 50 ul		100ea 74,600
P09-59-127	701910	for 100 ul		50ea 37,200
P09-59-129	702852	for 500 ul	PP ②	10ea 59,900
P09-59-131	702854	for 1ml		10ea 59,900
P09-59-133	702858	for 5ml		10ea 59,900
P09-59-135	702860	for 10ml		10ea 59,900

#### Seal

Cat. No.	Model	사용 용량	재질	PK Price ₩
P09-59-141	701920	10, 25 ul	PTFE ③	3ea 37,600
P09-59-143	701922	50 ul		3ea 37,600
P09-59-145	701924	100 ul		3ea 37,600
P09-59-147	702864	500 ul	PE ④	10ea 47,100
P09-59-149	702866	1 ml		10ea 47,100
P09-59-151	702870	5 ml		10ea 47,100
P09-59-153	702872	10 ml		10ea 47,100

#### Piston Rod with Seal

Cat. No.	Model	사용 용량	PK Price ₩
P09-59-161	701930	10 ul	3ea 65,800
P09-59-163	701932	25 ul	3ea 74,200
P09-59-165	701934	50 ul	3ea 74,200
P09-59-167	701936	100 ul ⑤	3ea 74,200

