



## "Brand" Pipette Tip 피펫 팁

- CE-marked according to EC-Directive 98/79 EC
- 재질 : PP (DiHEMA & Oleamide Free), Autoclavable at 121°C (20 min, 2bar)
- Free of DNA (<40fg) & RNase (<8.6 fg), Endotoxins (<1 pg) & ATP (<1 fg)  
: All Pipet Tips & Filter Tips
- Quick Volume Check를 위한 눈금을 가지고 있음 (Graduation)
- BIO-CERT® : Sterile, ISO11137, AAMI Guideline, USP29, Ph.Eur.Requirement  
- An SAL of 10<sup>-6</sup> is guaranteed
- Tip 사용 규격
  - 20 ~ 1000 ul Tip : Brand Pipette, Most of Pipette type of Gilson®, Thermo Fisher, Scientific Finnpiptette®, Eppendorf®, Biohit® / Sartorius®
  - 5 ml Tip : Brand Pipette, Thermo Fisher Scientific Finnpiptette®
  - 10 ml Tip : Brand Pipette, Eppendorf® & Gilson®
- Brand GmbH



### Bulk Tip

Cat. No.	Model	용량	눈금	길이 (mm)	PK Price ₩	CS Price ₩
P09-59-351	732002 (732022)	0.1~20 ul, Nano-cap	2, 20 ul	37	1000ea 41,100	10pk 384,100
P09-59-355	732004 (732024)	0.5 ~ 20 ul	2, 10 ul	46	1000ea 38,400	10pk 358,200
P09-59-357	732006 (732026)	1 ~ 50 ul	2.5, 10, 25, 50 ul	50	1000ea 28,000	10pk 261,200
P09-59-359	732008 (732008)	2 ~ 200 ul, Yellow Tip	20, 100 ul	50	1000ea 25,800	10pk 241,100
P09-59-363	732010 (732030)	5 ~ 300 ul	50, 100 ul	52	1000ea 34,200	10pk 319,600
P09-59-367	732012 (732032)	50 ~ 1000 ul, Blue Tip	250, 500 ul	70	1000ea 28,800	5pk 134,500
P09-59-375	702595	0.5 ~ 5 ml	X	160	200ea 43,700	
P09-59-379	702603	1 ~ 10 ml		156.5	100ea 40,800	

- P09-59-375 : Ø9.6 mm, P09-59-379 : Ø15 mm
- Subdivision Ref Model : 내부 소분 단위로 판매시 대응되는 "Brand" Model No. 임

### TipBox

- 재질 : PP(Polypropylene), Autoclavable at 121°C (20 min)
- Push-on Lid 포함, Stackable
- 5, 10 ml용은 원형임

#### Non-Sterile, 비멸균형 팁 박스 ①

Cat. No.	Model	용량	포장	PK Price ₩
P09-59-391	732202	0.1 ~ 20 ul	96 Tip/Box (5 Box/PK)	480ea 86,000
P09-59-395	732204	0.5 ~ 20 ul		480ea 80,500
P09-59-397	732206	1 ~ 50 ul		480ea 80,500
P09-59-399	732208	2 ~ 200 ul		480ea 76,900
P09-59-403	732210	5 ~ 300 ul		480ea 76,900
P09-59-407	732212	50 ~ 1000 ul		480ea 78,900
P09-59-415	702605	0.5 ~ 5ml	28 Tip/Box	140ea 385,300
P09-59-419	702608	1 ~ 10 ml	18 Tip/Box	90ea 498,000

#### BIO-CERT® Sterile, 멸균형 팁 박스 ②

Cat. No.	Model	용량	포장	PK Price ₩
P09-59-425	732222	0.1 ~ 20 ul	96 Tip/Box (10 Box/PK)	960ea 188,200
P09-59-429	732224	0.5 ~ 20 ul		960ea 175,300
P09-59-431	732226	1 ~ 50 ul		960ea 175,300
P09-59-433	732228	2 ~ 200 ul		960ea 161,400
P09-59-437	732230	5 ~ 300 ul		960ea 161,400
P09-59-441	732232	50 ~ 1000 ul		960ea 180,500

- 개별 멸균 포장



**Dual function hinged and push-on lid**  
Simply rotate 180° to change closure method

**Transparent window**  
For better all-round visual inspection

**Usable with a partially loaded multichannel pipette**  
The rimless tip-tray border enables problem-free loading of individual pipette tips with multichannel instruments

**No deflection during tip loading**  
The tip-tray is fabricated from especially rigid PP

**Colored tip-trays with side panel labeling**  
The contents of the box are always clearly visible

**Opening and closing with one hand**  
The clamping mechanism holds the tiptray securely in the box

다음 페이지에서 계속...