

## Precision Dispensing Tip, PD-Tip II 정밀형 분주기 팁

- With Tip Size Encoding
  - Handstep® Touch 에서 Tip 자가 인식을 위한 Encoding
- Batch Certificate : ISO8655 Requirement, CE-IVD-Compliant
- Eppendorf® Multipette® 4780/Repeater™, Gilson® Repetman®, Rainin AutoRep™ E,M,S와 같이 사용할 수 있음
- 재질 : PD Cylinder, PE-HD Piston (단, 0.1 ml : LCP)
- BIO-CERT®
  - Endotoxin-Free of DNA and RNase Free of ATP, Individually Wrapped
- Brand GmBH


### 비멸균형 Tip, Non-Sterile

Cat. No.	Model	Capacity	PK Price ₩
P09-59-932	705700	0.1 ml	100ea 154,000
P09-59-934	705702	0.5 ml	100ea 137,800
P09-59-936	705704	1.0 ml	100ea 137,800
P09-59-938	705706	1.25 ml	100ea 137,800
P09-59-940	705708	2.5 ml	100ea 137,800
P09-59-942	705710	5 ml	100ea 137,800
P09-59-944	705712	10 ml	100ea 137,800
P09-59-946	705714	12.5 ml	100ea 137,800
P09-59-948	705716	25 ml	50ea 77,600
P09-59-950	705718	50 ml	25ea 41,800
P09-59-952	705720	0.5 / 1 / 1.25 / 2.5 / 5 / 10 / 12.5 ml 각 20개	1Set 194,700

### 멸균형 Tip, BIO-CERT® Sterile

Cat. No.	Model	용량	PK Price ₩
P09-59-962	705730	0.1 ml	100ea 334,800
P09-59-964	705732	0.5 ml	100ea 295,500
P09-59-966	705734	1.0 ml	100ea 295,500
P09-59-968	705736	1.25 ml	100ea 295,500
P09-59-970	705738	2.5 ml	100ea 295,500
P09-59-972	705740	5 ml	100ea 295,500
P09-59-974	705742	10 ml	100ea 295,500
P09-59-976	705744	12.5 ml	100ea 295,500
P09-59-978	705746	25 ml	25ea 84,300
P09-59-980	705748	50 ml	25ea 87,200

- 개별포장 (Individually Wrapped)



### Tip Adapter for 25 ml & 50 ml Tip

Cat. No.	Model	Description	PK Price ₩
P09-59-991	702398	Non-Sterile Adapter	10ea 37,500
P09-59-993	702399	BIO-CERT® Sterile Adapter	1ea 8,100

### Volume Table Handystep® S with PD-Tips II

Set	Tip Size (ml)										Step
	0.1	0.5	1	1.25	2.5	5	10	12.5	25	50	
1	2	10	20	25	50	100	200	250	500	1000	49
1.5	3	15	30	37.5	75	150	300	375	750	1500	32
2	4	20	40	50	100	200	400	500	1000	2000	24
2.5	5	25	50	62.5	125	250	500	625	1250	2500	19
3	6	30	60	75	150	300	600	750	1500	3000	15
3.5	7	35	70	87.5	175	350	700	875	1750	3500	13
4	8	40	80	100	200	400	800	1000	2000	4000	11
4.5	9	45	90	112.5	225	450	900	1125	2250	4500	10
5	10	50	100	125	250	500	1000	1250	2500	5000	9
Volume (ul)											

## Repeating Pipette, Syringe Tip 보급형 주사기형 연속 분주기

Cat. No.	Model	분주 용량	사용 Tip	EA Price ₩
P09-92-233	H721008	10 ~ 5000 ul	0.5 ~ 50 ml	199,300

- 반복적인 피펫팅 작업용
- 농도가 있거나 점도가 있는 물질의 분주도 가능
- 한번 흡입으로 최대 48회 반복 분주 가능
- 7개의 PP Syringe Tip을 선택하여 사용
- Brand Dispensing Tip과 같이 사용하기에 알맞음 (상단 참조)
- DLAB

### Select the Dispensing Volume

Tips/Select	1	2	3	4	5
0.50 ml	10 ul	20 ul	30 ul	40 ul	50 ul
1.25 ml	25 ul	50 ul	75 ul	100 ul	125 ul
2.50 ml	50 ul	100 ul	150 ul	200 ul	250 ul
5.00 ml	100 ul	200 ul	300 ul	400 ul	500 ul
12.50 ml	250 ul	500 ul	750 ul	1000 ul	1250 ul
25.00 ml	500 ul	1000 ul	1500 ul	2000 ul	2500 ul
50.00 ml	1000 ul	2000 ul	3000 ul	4000 ul	5000 ul
Steps	48	23	15	11	8

