


Alumina Plate 알루미나 판 / 알루미나 플레이트


Cat. No.	Model	Size (w×d)	두께	EA Price ₩
P04-45-125	ALO0505	50×50 mm	5 mm	30,800
P04-45-132	ALO0510	50×50 mm	10 mm	41,000
P04-45-139	ALO1005	100×100 mm	5 mm	144,600
P04-45-146	ALO1010	100×100 mm	10 mm	162,700
P04-45-153	ALO1505	150×150 mm	5 mm	163,900
P04-45-160	ALO1510	150×150 mm	10 mm	204,900
P04-45-167	ALO2005	200×200 mm	5 mm	307,300
P04-45-174	ALO2010	200×200 mm	10 mm	469,900
P04-45-181	ALO3005	300×300 mm	5 mm	921,700
P04-45-188	ALO3010	300×300 mm	10 mm	994,000

- 재질 : Al₂O₃ (99.8%)
- 요청에 의해 규격의 사이즈 공급가능함
- LabCeramic®



Glass Plate 파이렉스 유리판 / 고품질 유리 플레이트

- Borosilicate Glass 3.3 고품질 유리판
- 광택 차수, 표면 평면도 우수함
- 열충격 및 급격한 온도 변화에 높은 저항력을 가짐
- 높은 온도에서 사용이 가능한 화학적 내구성을 지닌 특급 유리
- Max. Temp : 500 °C
- LK Labkorea



Circle Plate, 원형 유리판

Cat. No.	Model	Size	두께 (mm)	EA Price ₩
P04-132-658	1907-41	∅25 mm (1")	1.1	9,000
P04-132-662	1907-01	∅25 mm (1")	2.2	9,000
P04-132-666	1907-02	∅25 mm (1")	3.3	9,000
P04-132-670	1907-42	∅50 mm (2")	1.1	13,500
P04-132-674	1907-06	∅50 mm (2")	2.2	13,500
P04-132-678	1907-07	∅50 mm (2")	3.3	13,500
P04-132-682	1907-08	∅50 mm (2")	5	13,500
P04-132-686	1907-43	∅76 mm (3")	1.1	18,100
P04-132-690	1907-11	∅76 mm (3")	2.2	18,100
P04-132-694	1907-12	∅76 mm (3")	3.3	18,100
P04-132-698	1907-13	∅76 mm (3")	5	20,200
P04-132-702	1907-44	∅100 mm (4")	1.1	24,800
P04-132-706	1907-16	∅100 mm (4")	2.2	24,800
P04-132-710	1907-17	∅100 mm (4")	3.3	24,800
P04-132-714	1907-18	∅100 mm (4")	5	27,000
P04-132-718	1907-45	∅150 mm (6")	1.1	45,100
P04-132-722	1907-21	∅150 mm (6")	2.2	45,100
P04-132-726	1907-22	∅150 mm (6")	3.3	45,100
P04-132-730	1907-23	∅150 mm (6")	5	54,100
P04-132-734	1907-46	∅200 mm (8")	1.1	78,900
P04-132-738	1907-26	∅200 mm (8")	2.2	78,900
P04-132-742	1907-27	∅200 mm (8")	3.3	78,900
P04-132-746	1907-47	∅250 mm (10")	1.1	124,000
P04-132-750	1907-31	∅250 mm (10")	2.2	124,000
P04-132-754	1907-32	∅250 mm (10")	3.3	124,000
P04-132-758	1907-48	∅300 mm (12")	1.1	168,700
P04-132-762	1907-36	∅300 mm (12")	2.2	168,700
P04-132-766	1907-37	∅300 mm (12")	3.3	168,700

Square Plate, 사각 유리판

Cat. No.	Model	Size (mm)	두께 (mm)	EA Price ₩
P04-132-852	1905-35	25×25 (1")	1.1	9,000
P04-132-856	1905-55	25×25 (1")	2.2	9,000
P04-132-860	1905-66	25×25 (1")	3.3	9,000
P04-132-864	1905-36	50×50 (2")	1.1	13,500
P04-132-868	1905-01	50×50 (2")	2.2	13,500
P04-132-872	1905-02	50×50 (2")	3.3	13,500
P04-132-786	1905-37	76×76 (3")	1.1	18,100
P04-132-880	1905-06	76×76 (3")	2.2	18,100
P04-132-884	1905-07	76×76 (3")	3.3	18,100
P04-132-888	1905-38	100×100 (4")	1.1	24,800
P04-132-892	1905-11	100×100 (4")	2.2	24,800
P04-132-896	1905-12	100×100 (4")	3.3	24,800
P04-132-900	1905-13	100×100 (4")	5	27,000
P04-132-904	1905-39	150×150 (6")	1.1	45,100
P04-132-908	1905-16	150×150 (6")	2.2	45,100
P04-132-912	1905-17	150×150 (6")	3.3	45,100
P04-132-916	1905-18	150×150 (6")	5	54,100
P04-132-920	1905-40	200×200 (8")	1.1	78,900
P04-132-924	1905-21	200×200 (8")	2.2	78,900
P04-132-928	1905-22	200×200 (8")	3.3	78,900
P04-132-932	1905-41	250×250 (10")	1.1	124,000
P04-132-936	1905-26	250×250 (10")	2.2	124,000
P04-132-940	1905-27	250×250 (10")	3.3	124,000
P04-132-944	1905-42	300×300 (12")	1.1	168,700
P04-132-948	1905-31	300×300 (12")	2.2	168,700
P04-132-952	1905-32	300×300 (12")	3.3	168,700