


I Quartz Plate

석영 판

- Wavelength : 200 ~ 2500 nm
- 광학용 Quartz, 광택 차수 및 표면 평면도, 투과율 높음
- Max. Temp. : 1300°C, 상용 온도 : 1000°C
- YG Optical



Rectangular Quartz Plate, 사각 석영판

Cat. No.	Model	Size (w×d)	두께	EA Price ₩
P04-60-333	01UVC01	10×10 mm	1 T	6,100
P04-60-337	02UVC01	20×20 mm		9,200
P04-60-341	02UVC51	25×25 mm		9,800
P04-60-345	02UVC71	25×75 mm		18,200
P04-60-346	03UVC01	30×30 mm		13,000
P04-60-347	04UVC01	40×40 mm		16,200
P04-60-349	05UVC01	50×50 mm	1 T	21,300
P04-60-353	05UVC02		2 T	22,800
P04-60-357	05UVC03		3 T	24,400
P04-60-361	05UVC05		5 T	30,400
P04-60-365	05UVC09		10 T	41,000
P04-60-369	10UVC01	100×100 mm	1 T	65,000
P04-60-373	10UVC02		2 T	71,500
P04-60-377	10UVC03		3 T	78,000
P04-60-381	10UVC05		5 T	84,500
P04-60-385	10UVC09		10 T	129,100
P04-60-389	15UVC01	150×150 mm	1 T	123,500
P04-60-393	15UVC02		2 T	187,500
P04-60-397	15UVC03		3 T	136,500
P04-60-401	15UVC05		5 T	156,000
P04-60-405	15UVC09		10 T	260,000
P04-60-409	20UVC02	200×200 mm	2 T	240,500
P04-60-413	20UVC03		3 T	260,000
P04-60-417	20UVC05		5 T	286,000
P04-60-421	20UVC09		10 T	429,000
P04-60-425	25UVC02	250×250 mm	2 T	260,000
P04-60-429	25UVC03		3 T	299,000
P04-60-433	25UVC05		5 T	429,000
P04-60-437	25UVC09		10 T	624,000
P04-60-441	30UVC02	300×300 mm	2 T	364,000
P04-60-445	30UVC03		3 T	429,000
P04-60-449	30UVC05		5 T	598,000
P04-60-453	30UVC09		10 T	910,900

Circle Quartz Plate, 원형 석영판

Cat. No.	Model	Size	두께	EA Price ₩
P04-60-601	01UVR01	Ø10 mm	1 T	5,900
P04-60-603	01UVR02		2 T	6,500
P04-60-611	011UVR01	Ø15 mm	1 T	7,200
P04-60-613	011UVR02		2 T	7,900
P04-60-621	02UVR01	Ø20 mm	1 T	8,500
P04-60-623	02UVR02		2 T	9,100
P04-60-631	03UVR01	Ø30 mm	1 T	13,000
P04-60-633	03UVR02		2 T	14,300
P04-60-641	04UVR01	Ø40 mm	1 T	16,200
P04-60-643	04UVR02		2 T	18,200
P04-60-668	05UVR01	Ø50 mm	1 T	21,300
P04-60-672	05UVR02		2 T	22,800
P04-60-676	05UVR03		3 T	24,400
P04-60-680	05UVR05		5 T	30,400
P04-60-684	05UVR10		10 T	41,000
P04-60-688	10UVR01	Ø100 mm	1 T	65,000
P04-60-692	10UVR02		2 T	71,500
P04-60-696	10UVR03		3 T	78,000
P04-60-700	10UVR05		5 T	84,500
P04-60-704	10UVR10		10 T	129,100
P04-60-712	15UVR02	Ø150 mm	2 T	130,000
P04-60-716	15UVR03		3 T	136,500
P04-60-720	15UVR05		5 T	156,000
P04-60-724	15UVR10	Ø200 mm	10 T	325,000
P04-60-728	20UVR02		2 T	240,500
P04-60-732	20UVR03		3 T	260,000
P04-60-736	20UVR05	Ø250 mm	5 T	286,000
P04-60-740	20UVR10		10 T	429,000
P04-60-744	25UVR02	Ø250 mm	2 T	260,400
P04-60-748	25UVR03		3 T	296,800
P04-60-752	25UVR05		5 T	363,800
P04-60-756	25UVR10		10 T	520,700

I Silicon Plate

실리콘 플레이트, 300×300mm

Cat. No.	Model	Size (w×d)	두께	색상	EA Price ₩
P04-140-105	LS0005	300×300 mm	0.5 T	Natural	5,300
P04-140-111	LS0010		1.0 T		5,200
P04-140-139	LS0015		1.5 T		8,700
P04-140-141	LS0200		2.0 T		12,300
P04-140-171	LS0300		3.0 T		18,500
P04-140-201	LS0500		5.0 T		29,300
P04-140-221	LS0700		7.0 T		32,900
P04-140-239	LS1000		10.0 T		52,100

- 재질 : 실리콘 고무
- 내열, 내한, 내유성, 내약품성이 뛰어남
- 내열 온도 : 200°C, 내한 온도 : -50°C
- 경도 : 50°
- Tommyheco

