


Large Silicon Plate 대형 실리콘 플레이트

Cat. No.	Model	Size	두께	EA Price ₩
P04-140-525	LS1000	1 × 1 m	1 T	38,600
P04-140-527	LS1001		1.5 T	62,700
P04-140-529	LS1011		2 T	77,200
P04-140-532	LS1002		3 T	115,700
P04-140-539	LS1003		5 T	195,200


- Shore (경도) : 50°, Natural Color
- Tommyheco



PTFE Sheet Plate PTFE 테프론 시트 플레이트

Cat. No.	Model	Size	두께	EA Price ₩
P04-140-777	F8034-01	300 mm × 10 m	50 um	300,200
P04-140-779	F8034-02		80 um	190,900
P04-140-782	LP13100		100 um	125,400
P04-140-806	LP13200		200 um	120,300
P04-140-812	LP13300	300 mm × 1 m	300 um	25,100
P04-140-816	LP13500		500 um	41,900
P04-140-822	LP13800		800 um	54,900
P04-140-828	LP14100		1.0 mm	77,000

- 재질 : PTFE • 내열 온도 : 260°C • 내약품, 저마찰계수, 비점착성 등이 우수함
- Tommyheco



PTFE Sheet Plate PTFE 테프론 비접착식 시트 플레이트

Cat. No.	Model	Size	두께	EA Price ₩	PK Price ₩
P04-460-101	S1803-02	300 mm × 1 m	0.05 mm	-	3ea 66,900
P04-460-104	S1803-04		0.12 mm	-	3ea 80,000
P04-460-107	S1803-06		0.25 mm	45,200	-
P04-460-110	S1803-08		0.5 mm	79,300	-
P04-460-113	S1803-10		0.75 mm	112,200	-
P04-460-116	S1803-12		1.0 mm	126,900	-
P04-460-119	S1803-14		1.5 mm	196,600	-
P04-460-122	S1803-21	600 mm × 1 m	0.05 mm	-	3ea 91,400
P04-460-125	S1803-23		0.12 mm	38,200	-
P04-460-128	S1803-25		0.25 mm	78,900	-
P04-460-131	S1803-27		0.5 mm	144,300	-
P04-460-134	S1803-29		0.75 mm	213,700	-
P04-460-137	S1803-31		1.0 mm	241,600	-
P04-460-140	S1803-33		1.5 mm	379,300	-

- 재질 : PTFE, -200 °C to 250°C 내화학 및 내열 • 색상 : White
- Table Pad, Lining Drawer, Slideway, Insulation 용도에 적합
- BOLA



FEP Sheet Plate FEP 테프론 비접착식 시트 플레이트

Cat. No.	Model	용량	두께	EA Price ₩
P04-460-010	S1833-04	150 mm × 1 m	0.025 mm	79,300
P04-460-013	S1833-08		0.05 mm	72,300
P04-460-016	S1833-16		0.25 mm	345,100
P04-460-019	S1833-34	300 mm × 1 m	0.025 mm	90,100
P04-460-022	S1833-38		0.05 mm	111,400
P04-460-025	S1833-46		0.25 mm	670,700

- 재질 : FEP, -200 °C to 205°C 내화학 및 내열
- Transparent, Gastight and Non-porous rolls
- Table Pad, Lining Drawer, Slideway, Insulation 용도에 적합
- BOLA

