

BRANDplate® Cell Culture Plate & Insert 셀 컬처 플레이트와 인서트

- Manual & Automated Cell & Tissue Culture용
- 24 Well Standard Plate & 6 Well Special Plate (PS재질, Sterile)
- pureGrade™ S (Untreated) & cellGrade™ plus (Hydrophilic Cell Culture)
- Side Feed Port in the Well for Changing Media with Minimal Disruption
- Perfect for Use with Pipetting Robot
- Especially for Tissue Culture : Inserts with Inlet Channel
- Manufactured according to Requirements in ANSI/SLAS Standard 1 & 4
- Non-cytotoxic according to ISO 10 993-5, Endotoxins-free (< 0.01 EU/mL)
- Free from DNA, DNase & RNase
- Sterile according to ISO 11137 & AAMI Guidelines, SAL 10⁻⁶
- Brand GmbH

BRANDplates® Microplate, Sterile, with Lid & Condensation Ring, 개별포장

Cat. No.	Model	Description	Surface	PK Price ₩
P04-59-761	782880	24-well Standard	pureGrade™ S (Un-treated)	10ea 213,700
P04-59-763	782881	6-well		10ea 213,700
P04-59-765	782890	24-well Standard	cellGrade™ plus	10ea 224,300
P04-59-767	782891	6-well		10ea 224,300

- 24-well Standard Plate : 24 Individually Fillable Well, 4 Insert Strip 또는 개별 Insert 와 사용 가능
- 6-well Special Plate : 4 Well이 하나로 연결된 형태, 4 Insert Strip 과 사용하기에 알맞음

BRANDplates® Insert Strip

- PS, cellGrade™ plus Surface, Sterile, Strips of 4 Insert (Divisible), 개별포장

Cat. No.	Model	Type	Height	Membrane	Pore	PK Price ₩
P04-59-771	782800	Smooth-walled	13 mm	PC	0.4 um	12ea 488,000
P04-59-772	782820				1 um	12ea 488,000
P04-59-773	782840				3 um	12ea 488,000
P04-59-774	782860				8 um	12ea 488,000
P04-59-795	782900				9 mm	0.4 um
P04-59-775	782810	13 mm	PET	0.4 um	12ea 488,000	
P04-59-776	782870			8 um	12ea 488,000	
P04-59-777	782801	with Inlet Opening System	13 mm	PC	0.4 um	12ea 585,300
P04-59-778	782821				1 um	12ea 585,300
P04-59-779	782841				3 um	12ea 585,300
P04-59-780	782861		8 um	12ea 585,300		
P04-59-796	782901		9 mm	0.4 um	12ea 646,900	
P04-59-781	782811		13 mm	PET	0.4 um	12ea 585,300
P04-59-782	782871				8 um	12ea 585,300

- Hanging Insert : 4mm Distance (from Well Bottom to Insert Membrane)
- Inlet Opening System : Automated in Vitro Preparation of Human Tissues
 - Rapid, Consistent Changing of Media, from Submersion Culture to Air-lift Culture
 - Adjustment of Medium Level without Damage to Skin Model


BRANDplates® Single Insert

- PS, cellGrade™ plus Surface, Sterile, Individual Insert

Cat. No.	Model	Type	Height	Membrane	Pore	PK Price ₩
P04-59-870	782806	Smooth-walled	13 mm	PC	0.4 um	48ea 558,200
P04-59-873	782826				1 um	48ea 563,700
P04-59-876	782846				3 um	48ea 563,700
P04-59-879	782866			PET	8 um	48ea 563,700
P04-59-882	782816				0.4 um	48ea 563,700
P04-59-885	782876				8 um	48ea 563,700



24-well standard plate



6-well special plate



Insert strips, smooth-walled



Insert strips with special inlet channels (the Inlet Opening System*)

Pore size	Membrane / Pore size Application Example
0.4 um	Co-culture, Transport Study, Secretion Study, Cell Polarity Study, etc.
1.0 um	Co-culture, Transport Study, Secretion Study, etc.
3.0 & 8.0 um	Migration Studies, Chemotaxis Study, Metastasis Experiment, etc.