


96-Well Microplate, Pure Grade™ 96 웰 플레이트, Non-Treated, Non-Sterile, ELISA Assay

Cat. No.	Model	재질	Type	색상	Bottom Type	용량	PK Price ₩		
P04-59-111	781600	PS	Standard	Transparent	U	330 ul	100ea	180,400	
P04-59-114	781601				V	360 ul	100ea	180,400	
P04-59-117	781602				F	100ea	180,400		
P04-59-120	781603				C	350 ul	100ea	180,400	
P04-59-123	781604			White	U	330 ul	100ea	483,600	
P04-59-126	781605				F	350 ul	100ea	483,600	
P04-59-129	781607				Black	U	330 ul	100ea	483,600
P04-59-132	781608					F	350 ul	100ea	483,600
P04-59-135	781610		Transp. Bottom	White	F	330 ul	100ea	1,289,000	
P04-59-138	781611						Black	100ea	1,289,000
P04-59-141	781614		UV-Transparent	Transparent	F	350 ul	50ea	1,207,000	


- 비열균, 표면 비처리, 대부분의 실험에 적용 가능한 범용 마이크로 플레이트
- 균일 분석용(단백질, 핵산등의 분석)이나 HTS를 위한 Sample 저장용으로 특화됨
- 화이트, 블랙 플레이트의 경우 회색 양각된 Alphanumeric 코드로 명확한 구별이 가능하게 제작됨
- Free from Endotoxins, DNase, DNA, RNase, Non-cytotoxic
- UV Polymer 재질의 마이크로 플레이트는 낮은 파장에서도 높은 투과도를 보여 단백질, 핵산 분석에 유리하며 화학적인 내성이뛰어남 (하단 자료 참조)
- Brand GmBH

Resistance of Frequently Used Plastics

Media	UV-Polymer	PS	PP
Acetic Acid, 100%	+	-	+
Acetone	+	-	+
Ammonia	+	+	+
Benzene(Benzol)	-	-	+
Chloroform	-	-	-
Dioxane	+	-	+
Ethyl Acetate	+	-	+
Hexane	-	-	+
Hydrochloric Acid, 36%	+	+	+
Hydrofluoric Acid, 10%	+	+	+
Isopropanol	+	+	+
Methanol	+	+	+



Transmittion of PS and UV Polymer




- UV Polymer : Slightly Hydrophilized, Bottom UV-transparent Film, 25 um Wall Thickness, Skirt PS, Transmission at 240 nm 80%

96-Well Microplate with Lid, PureGrade™ S 96 웰 플레이트, Non-Treated, Sterile, ELISA Assay

Cat. No.	Model	재질	Type	색상	Bottom Type	용량	PK Price ₩				
P04-59-171	781660	PS	Standard	Transparent	U	330 ul	50ea	230,400			
P04-59-174	781661				V	360 ul	50ea	230,400			
P04-59-177	781662				F	50ea	230,400				
P04-59-180	781663				C	50ea	230,400				
P04-59-183	781665			White	F	350 ul	50ea	413,700			
P04-59-186	781668						Black	50ea	413,700		
P04-59-189	781670						Transp. Bottom	White	330 ul	50ea	1,100,200
P04-59-192	781671									Black	50ea

- 열균(β-radiation), 표면 비처리
- Bacteriological Assays에 특화됨
- 화이트, 블랙 플레이트의 경우 회색 양각된 Alphanumeric 코드로 명확한 구별이 가능하게 제작됨
- Free from Endotoxins, DNase, DNA, RNase, Non-cytotoxic
- Brand GmBH



"Brand" 사의 다양한 96-well Microplate는 www.labmart.co.kr 참조