

## Digital Manometer / Pressure / Vacuum Gauge 디지털 차압계, Traceable®

Cat. No.	Model	압력범위 (PSI)	진공범위 (mbar)	EA Price ₩
V03-147-010	3460	-5.000 ~ 5.000	-345.0 ~ 345.0	464,500
V03-147-020	3461	-15.00 ~ 15.00	-1034 ~ 1034	464,500
V03-147-030	3462	-30.00 ~ 30.00	-2068 ~ 2068	464,500
V03-147-040	3463	-100.0 ~ 100.0	-6895 ~ 6895	464,500
V03-147-050	3466	-5.000 ~ 5.000	-345.0 ~ 345.0	1,043,600
V03-147-060	3469	-100.0 ~ 100.0	-6895 ~ 6895	1,043,600

- 11개의 압력, 진공 단위를 선택해서 표시 가능
- 메모리 기능으로 최대 & 최소 표시 확인 가능 및 Hold Key 표시값 고정 가능
- Display 백라이트로 어두운 곳에서 관찰 가능
- 프린터 가능제품(3466,3469)은 샘플번호, 날짜, 시간, 현재값 표시하고 1초~2시간 간격 인쇄 설정 가능
- 2가지 튜브(직경 4mm or 1/6" x 길이 483mm) , 9V 건전지, 케이스 포함
- Count Timer : up to 99 hour 59 min 59 sec
- 크기 : w70 x d29 x h181mm
- 국내 소급 가능한 Traceable® 성적서 제공 (ILAC)
- Cole-Parmer / Traceable®



3460



3466

Cat. No.	V03-147-010	V03-147-020	V03-147-030	V03-147-040	V03-147-050	V03-147-060
Model	3460	3461	3462	3463	3466	3469
최대허용 압력(PSI)	20	30	60	150	20	150
PSI Range	-5.000 to 5.000	-15.00 to 15.00	-30.00 to 30.00	-100.0 to 100.0	-5.000 to 5.000	-100.0 to 100.0
PSI Resolution	0.003	0.01	0.02	0.1	0.003	0.1
mbar Range	-345.0 to 345.0	-1034 to 1034	-2068 to 2068	-6895 to 6895	-345.0 to 345.0	-6895 to 6895
mbar Resolution	0.2	1	2	4	0.2	4
bar Range	-0.3450 to 0.3450	-1.034 to 1.034	-2.068 to 2.068	-6.895 to 6.895	-0.3450 to 0.3450	-6.895 to 6.895
bar Resolution	0.001	0.001	0.002	0.004	0.001	0.004
cm/H2O Range	-351.0 to 351.0	-1055 to 1055	-2109 to 2109	-7031 to 7031	-351.0 to 351.0	-7031 to 7031
cm/H2O Resolution	0.2	1	2	4	0.2	4
inch/H2O Range	-138.0 to 138.0	-415.0 to 415.0	-830.0 to 830.0	-2768 to 2768	-138.0 to 138.0	-2768 to 2768
inch/H2O Resolution	0.1	0.3	0.5	2	0.1	2
feet/H2O Range	-12.00 to 12.00	-35.00 to 35.00	-69.00 to 69.00	-231.0 to 231.0	-12.00 to 12.00	-231.0 to 231.0
feet/H2O Resolution	0.01	0.02	0.04	0.2	0.01	0.2
mm/Hg Range	-259.0 to 259.0	-776.0 to 776.0	-1551 to 1551	-5171 to 5171	-259.0 to 259.0	-5171 to 5171
mm/Hg Resolution	0.2	0.5	1	3	0.2	3
inch/Hg Range	-10.00 to 10.00	-31.00 to 31.00	-61.00 to 61.00	-204.00 to 204.00	-10.00 to 10.00	-204.00 to 204.00
inch/Hg Resolution	0.01	0.03	0.1	1	0.01	1
kg/cm2 Range	-0.350 to 0.350	-1.060 to 1.060	-2.110 to 2.110	-7.030 to 7.030	-0.350 to 0.350	-7.030 to 7.030
kg/cm2 Resolution	0.001	0.001	0.002	0.004	0.001	0.004
oz/inch2 Range	-80.00 to 80.00	-240.0 to 240.0	-480.0 to 480.0	-1600 to 1600	-80.00 to 80.00	-1600 to 1600
oz/inch2 Resolution	0.05	0.2	0.3	1	0.05	1
kPa Range	-34.50 to 34.50	-103.0 to 103.0	-207.0 to 207.0	-690.0 to 690.0	-34.50 to 34.50	-690.0 to 690.0
kPa Resolution	0.02	0.1	0.2	0.4	0.02	0.4

### Optional Accessory

Cat. No.	Model	Description	EA Price ₩
V03-147-071	7337	Replacement Paper rolls (3466 & 3469 only)	13,400

- Analog Manometer는 1203 Page 참조
- Digital Vacuum Gauge는 1201 ~ 1203 Page 참조

