


## Manual Liquid Pump 수동식 액체 이송 펌프 / 자바라

Cat. No.	Model	토출량	본체 (mm)	Arm	EA Price ₩
P06-206-034	LP-M01	7L / min	소, Ø14xh400	550 mm	1,800
P06-206-037	LP-M11	12L / min	중, Ø25xh610	660 mm	4,700
P06-206-040	LP-M02	17.5L / min	대(드럼용), Ø22xh870	1200 mm	12,200


- 재질 : PE (Polyethylene)
- M11은 원터치 스톱 누름봉으로 멈춤이 쉬움
- 드럼펌프(M02) 재질 : Body-PE, Pipe-PVC
- LK Labkorea



## Manual Acid-Liquid Pump 수동식 내산 액체 이송 펌프 / 자바라

Cat. No.	Model	토출량	본체 (mm)	Arm	EA Price ₩
P06-206-051	LP-M01S	7L / min	소, Ø15xh530	540 mm	2,800
P06-206-054	LP-M11S	12L / min	중, Ø25xh610	660 mm	5,300
P06-206-057	LP-M02S	21L / min	대(드럼용), Ø25.5xh1150	1200 mm	13,300

- 재질 : PE (Polyethylene)
- M11S는 원터치 스톱 누름봉으로 멈춤이 쉬움
- 등유, 가솔린, 디젤, 톨루엔, 시너(희석제), 질산 25%, 염산 35%, 황산 35%, 초산 35%, 물 등
- LK Labkorea



## Drum Pump 대용량 레버식 펌프

Cat. No.	Model	토출량	본체 (mm)	Arm	EA Price ₩
P06-206-061	LP-M53	530ml/stroke	Ø25 x H690-980	2m	39,600


- 용도 : 고점도 오일, 등유, 경유, 유압유(crt1000이하) 등
- 스트로크 당 530ml 토출 가능
- 파이프 높이 조절이 가능하여 드럼통에서 사용 가능
- LK Labkorea



## Small Liquid Pump 소형 액체 이송 펌프

Cat. No.	Model	용량	Tube Size	EA Price ₩
P06-20-305	9001210	4 L/min	Ø10×h450 mm	127,300
P06-20-315	9001112	14 L/min	Ø12×h700 mm	174,200
P06-20-325	9001115	20 L/min	Ø15×h700 mm	182,900

- 재질 : PP Delivery Tube, PVC Stopper
- 사용 용기 구성 : Ø30 ~ 70 mm (단, 9001210은 Ø25 ~ 36 mm입)
- 다양한 사이즈의 Bottle, Flask용 Handy Pump임
- Wenk Labtech



## PP Foot Pump PP 풋 펌프

- 재질 : PP, -10 ~ 125°C 내열
- PVC Barrel : Ø30 ~ Ø70 mm용
- Universal Bungs for Container Openings of Ø30 ~ Ø70 mm
- 내화학성이 강한 PP 재질로 Acid, Alkali, Aggressive Detergent 사용에 적합
- Wenk Labtec

with PP Delivery Tube ①

Cat. No.	Model	Flow Rate	Immersion Tube	EA Price ₩
P06-100-205	9001018	26 L/min	Ø18 mm x L700	221,000

with PVC Delivery Tube ②

Cat. No.	Model	Flow Rate	Immersion Tube	EA Price ₩
P06-100-255	9001133	12 L/min	Ø12 mm x L700	221,000

