

Pump Head for Peristaltic Pump 정량 펌프용 헤드, Masterflex®

• Masterflex®

Standard Pump Head ①

Cat. No.	Model	Tubing Size	EA Price ₩
P06-55-832	7013-20	LS13	384,600
P06-55-836	7014-20	LS14	384,600
P06-55-840	7016-20	LS16	384,600
P06-55-844	7017-20	LS17	384,600
P06-55-848	7018-20	LS18	384,600
P06-55-852	7015-20	LS15	384,600
P06-55-856	7024-20	LS24	384,600
P06-55-860	7035-20	LS35	384,600
P06-55-864	7036-30	LS36	513,800

- 재질 : PC Housing & CRS (Cold-rolled Steel) Rotor
- Silicon Tubing 기본 포함
- Driver에 따라 4ea까지 연결하여 사용가능 (Mounting Hardware는 별도임)

Pump Head, Easy-Load® III ②

Cat. No.	Model	Tubing Size	Rotor	EA Price ₩
P06-55-930	77800-50	LS13, LS14, LS16,	CRS	612,400
P06-55-934	77800-60	LS25, LS17, LS18	SS	781,900
P06-55-938	77800-52	LS15, LS24, LS35,	CRS	612,400
P06-55-942	77800-62	LS36	SS	781,900

- 재질 : PP Housing & 3ea (CRS or SS (Stainless Steel)) Rotor
- Mounting Hardware 없이 Head 끼리 연결 가능
- Mounting Plate/Adapter Kit와 Silicon Tubing 기본 포함
- Driver에 따라 4ea까지 연결하여 사용가능

Pump Head, Easy-Load® II ③

Cat. No.	Model	Tubing Size	Rotor	EA Price ₩
P06-55-914	77200-50	LS13, LS14, LS16,	CRS	664,600
P06-55-917	77200-60	LS25, LS17, LS18	SS	838,800
P06-55-920	77200-52	LS15, LS24, LS35,	CRS	664,600
P06-55-923	77200-62	LS36	SS	838,800
P06-55-925	77202-50	LS13, LS14, LS16,	CRS	921,700
P06-55-927	77202-60	LS25(2-Channel)	SS	914,300

- 재질 : PPS (Polyphenylene Sulfide) Housing & 4ea (CRS Rotor or SS) Rotor
- Drive 연결용 Mounting Hardware (Single Channel)와 Silicon Tubing 기본 포함
- Driver에 따라 4ea까지 연결하여 사용가능 (Mounting Hardware는 별도임)

Pump Head, Easy-Load® I ④

Cat. No.	Model	Tubing size	Rotor	EA Price ₩
P06-55-991	07514-10	LS13, LS14, LS16, LS25, LS17, LS18	SS	719,800
P06-55-292	07514-12	LS15, LS24	SS	719,800

- Driver에 따라 4ea까지 연결하여 사용가능 (Mounting Hardware는 별도임)

Mounting Hardware ⑤

Cat. No.	Model	Pump Head 수	적용 Pump Head	EA Price ₩
P06-55-950	7013-04	1 ea	Standard(LS36 제외) Easy-Load® I	77,000
P06-55-953	7013-05	2 ea		82,000
P06-55-956	7013-08	3 ea		95,900
P06-55-959	7013-09	4 ea		106,000
P06-55-962	7036-01	1 ea	Standard LS36 용	123,500
P06-55-965	7036-02	2 ea		183,700
P06-55-968	77200-01	1 ea	Standard Easy-Load® II	73,800
P06-55-971	77200-02	2 ea		84,500
P06-55-974	77200-03	3 ea		106,000
P06-55-977	77200-04	4 ea		122,500

- Pump Head를 다단으로 장착할때 쓰임



Standard Pump Head



Easy-Load® III



Easy-Load® II



Easy-Load® I

