

## Peristaltic Pump Tubing 정량 펌프용 튜빙, Masterflex®



Platinum Silicon



PharMed BPT



Chem-Durance® Bio



Viton

재질	Peroxide-cured Silicon / Platinum Silicon	PharMed BPT	Chem-Durance® Bio	Viton	Tygon® E-LFL	Tygon® E-Lab
내열온도	-50 ~ 230°C	-51 ~ 132°C	-60 ~ 74°C	-32 ~ 205°C	-46 ~ 74°C	-46 ~ 74°C
저항성	Acid	Not Recommended	Good	Excellent	Good	Good
	Alkalies	Not Recommended	Good	Excellent	Excellent	Good
	Organic Solvent	Not Recommended	Not Recommended	Good	Variable-Test before Using	Not Recommended
	Pressure	Fair	Good	Excellent	Good	Good
	Vacuum	Good	Good	Excellent	Good	Good
	Viscous Fluids	Fair	Excellent	Excellent	Good	Excellent
	Sterile Fluids	Excellent	Good	Excellent	Fair	Good

### Precision Pump Tubing

재질		LS13	LS14	LS16	LS25	LS17	LS18
Platinum Silicon (7.6m, 25ft)	Cat. No.	T23-55-011	T23-55-015	T23-55-019	T23-55-023	T23-55-031	T23-55-039
	Model	96410-13	96410-14	96410-16	96410-25	96410-17	96410-18
	EA Price ₩	122,400	150,000	144,600	189,500	200,500	234,500
PharMed-BPT (7.6m, 25ft)	Cat. No.	T23-55-111	T23-55-115	T23-55-119	T23-55-123	T23-55-131	T23-55-139
	Model	6508-13	6508-14	6508-16	6508-25	6508-17	6508-18
	EA Price ₩	186,900	241,300	241,300	250,900	302,800	309,600
Chem-Durance Bio® (15.2m, 50ft)	Cat. No.	T23-55-212	T23-55-216	T23-55-220	T23-55-224	T23-55-232	T23-55-240
	Model	6442-13	6442-14	6442-16	6442-25	6442-17	6442-18
	EA Price ₩	334,000	334,000	454,000	526,300	625,800	749,800
Viton (7.6m, 25ft)	Cat. No.	T23-55-311	T23-55-315	T23-55-319	T23-55-323	T23-55-331	T23-55-339
	Model	96412-13	96412-14	96412-16	96412-25	96412-17	96412-18
	EA Price ₩	227,700	244,000	305,500	398,200	508,500	586,300
Tygon® E-LFL (7.6m, 25ft)	Cat. No.	T23-55-421	T23-55-424	T23-55-427	T23-55-430	T23-55-436	T23-55-433
	Model	6440-13	6440-14	6440-16	6440-25	6440-17	6440-18
	EA Price ₩	159,600	195,000	212,800	225,000	259,100	291,800
Tygon® E-Lab (15.2m, 50ft)	Cat. No.	T23-55-451	T23-55-454	T23-55-457	T23-55-460	T23-55-466	T23-55-463
	Model	6509-13	6509-14	6509-16	6509-25	6509-17	6509-18
	EA Price ₩	90,900	105,700	154,200	211,300	248,200	332,800

Specification	ID	0.8 mm	1.6 mm	3.1 mm	4.8 mm	6.4 mm	7.9 mm
	Hose Barb	1.6 mm	1.6 mm	3.2 mm	4.8 mm	6.4 mm	9.5 mm
	최대 압력	40 psig			35 psig	20 psig	15 psig

### High Precision Pump Tubing

재질		LS15	LS24	LS35	LS36
Peroxide-cured Silicon (7.6m, 25ft)	Cat. No.	T23-55-501	T23-55-504	T23-55-507	T23-55-510
	Model	96400-15	96400-24	96400-35	96400-36
	EA Price ₩	192,300	244,000	306,800	335,500
Platinum Silicon (7.6m, 25ft)	Cat. No.	T23-55-527	T23-55-535	T23-55-543	T23-55-551
	Model	96410-15	96410-24	96410-35	96410-36
	EA Price ₩	196,400	267,200	317,600	362,600
PharMed-BPT (7.6m, 25ft)	Cat. No.	T23-55-627	T23-55-635	T23-55-643	T23-55-651
	Model	6508-15	6508-24	6508-35	6508-36
	EA Price ₩	241,300	334,000	368,200	381,800
Chem-Durance Bio® (15.2m, 50ft)	Cat. No.	T23-55-728	T23-55-736	T23-55-744	T23-55-752
	Model	6442-15	6442-24	6442-35	6442-36
	EA Price ₩	614,900	846,600	1,044,300	1,041,500
Viton (7.6m, 25ft)	Cat. No.	T23-55-827	T23-55-835	T23-55-843	T23-55-851
	Model	96412-15	96412-24	96412-35	96412-36
	EA Price ₩	642,200	768,900	963,800	1,025,100

Specification	ID	4.8 mm	6.4 mm	7.9 mm	9.7 mm
	Hose Barb	4.8 mm	6.4 mm	9.5 mm	9.5 mm
	최대 압력	40 psig		35 psig	20 psig

\* www.labmart.co.kr에서 좀 더 많은 제품을 확인할 수 있습니다.