

Industrial Pump Tubing 산업용 펌프 튜빙, Masterflex®



구분	Description
Silicon (Platinum-cured) ①	- Complies with USP Class V, EP, FDA, REACH, RoHS, ADCF - Fungus-Resistant; Phthalate- and Latex-free - Accommodates a wide Temperature Range
C-Flex® ②	- Complies with USP Class VI, EP, REACH, RoHS - Heat Sealable, Weldable, Economical - Biocompatible
Tygon® E-LFL ③	- Complies with PUS Class VI, EP, FDA, EU Food, REACH, RoHS, ABCF - Longest Life of all Tygon Formulations - Best Chemical Resistance
Tygon Chemical ④	- Complies with FDA 21 CFR 175.300, REACH, RoHS, ABCF - Best Chemical Resistance of any Tygon Formulation - Plasticizer-free and Low Gas Permeability
Viton ⑤	- Complies with FDA, ABCF - Perfect where FDA and Chemical Compatibility are needed - Resists Corrosives, Solvents, and Oils at Elevated Temperatures
Puri-Flex™ ⑥	- Complies with USP Class VI, EP, REACH, RoHS, ADCF - Heat Sealable, Weldable, Economical - Long Pump Life when compared to Silicone or many other TPE Tubings
Norprene ⑦	- Complies with NSF, REACH, RoHS - Up to 10,000 hours of Tubing Life - Best Choice for Pressure / Vacuum Applications

Precision Pump Tubing

재질		IP26	IP73	IP82
Silicon (Platinum-cured)	Cat. No.	P06-55-511	P06-55-515	P06-55-519
	Model	96410-26	96410-73	96410-82
	EA Price ₩	355,900	394,000	583,500
C-Flex®	Cat. No.	P06-55-771	P06-55-774	P06-55-777
	Model	06424-26	06424-73	06424-82
	EA Price ₩	263,200	246,900	310,900
Puri-Flex™	Cat. No.	P06-55-659	P06-55-663	P06-55-667
	Model	96419-26	96419-73	96419-82
	EA Price ₩	231,800	237,300	299,600
Tygon® E-LFL	Cat. No.	P06-55-535	P06-55-539	P06-55-543
	Model	06440-26	06440-73	06440-82
	EA Price ₩	449,900	432,200	563,100
Tygon Chemical	Cat. No.	P06-55-547	P06-55-551	P06-55-555
	Model	06475-26	06475-73	06475-82
	EA Price ₩	489,500	718,400	732,000
Norprene	Cat. No.	P06-55-671	P06-55-675	P06-55-679
	Model	06404-26	06404-73	06404-82
	EA Price ₩	331,300	396,800	522,200
Viton	Cat. No.	P06-55-559	P06-55-563	-
	Model	96412-26	96412-73	-
	EA Price ₩	1,183,200	1,237,800	-
Specification	ID			
	Hose Barb	6.4 mm	9.5 mm	12.7 mm
	Flow Range	0.01~4 LPM	0.01~8 LPM	0.02~13 LPM

High-Performance Pump Tubing

재질		IP70	IP88	IP89
Silicon (Platinum-cured)	Cat. No.	P06-55-611	P06-55-615	P06-55-619
	Model	96510-70	96510-88	96510-89
	EA Price ₩	564,400	635,300	833,000
Puri-Flex™	Cat. No.	P06-55-623	P06-55-627	P06-55-631
	Model	96419-70	96419-88	96419-89
	EA Price ₩	197,800	211,300	257,700
Tygon® E-LFL	Cat. No.	P06-55-635	P06-55-639	P06-55-643
	Model	06440-70	06440-88	06440-89
	EA Price ₩	793,500	980,200	1,025,100
Norprene	Cat. No.	P06-55-647	P06-55-651	P06-55-655
	Model	06404-70	06404-88	06404-89
	EA Price ₩	407,700	573,900	642,200
Specification	ID			
	Hose Barb	9.5 mm	12.7 mm	15.88 mm
	Flow Range	0.01~8 LPM	0.02~17 LPM	0.03~19 LPM