

IKA Chemical Resistant Vacuum Pump 내화학용 진공 펌프, MVP 10 Basic

Cat. No.	Model	유량	최대진공도	EA Price
P06-05-001	MVP10 Basic	28.3 L/min	7 mbar	3,360,000

- For Dry, Oil-free application
- 내화학성이 우수한 Diaphragm Pump
- Gas와 접촉하는 부분들은 PTFE 화합물로 구성 • Compact한 디자인 및 긴 수명
- Protection class according to DIN EN 60529 (IP54) • IKA

Model		MVP10 Basic
Pumping Speed		28.3 L/min
Ultimate Pressure	w/Gas Ballast	15 mbar
	w/o Gas Ballast	7 mbar
Inlet/Outlet		NW 10
Diaphragm Material		PTFE
Valve Material		EFPM
Noise Level		65 dB(A)
Size/Weight		w245×d320×h155 mm/11.5 kg
Environmental Conditions		Temp. 10-40°C / Humidity 80%
Power Input		180W
Power		230 V, 50/60 Hz

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
P06-05-011	20005194	VSS 1 Vacuum Safety Set	1,780,000
P06-05-007	20004850	VCV 1 Manual Vacuum Control Valve (Stand 별도)	1,170,000
P06-05-009	20004893	VCV 2 Digital Vacuum Control Valve (Stand 별도)	2,020,000
P06-05-017	20013886	VH.SI.8 Vacuum Silicone hose, Ø8 mm	80,000

IKA Chemical Resistant Diaphragm Vacuum Pump 다이아프램 진공 펌프, VACSTAR

Cat. No.	Model	유량	최대진공도	EA Price
P06-05-003	VACSTAR Digital	22 L/min	2 mbar	4,970,000
P06-05-004	VACSTAR Control			7,820,000

- 4개의 Chamber Membrane으로 이루어진 Vacuum Pump
- 구성 : Digital Type - 본품, Vacuum Hose(1m), Silencer, USB A-B Cable, Connecting Cable
Control Type - 본품, WiCo Control, Vacuum Hose(1m), Wico Bracket, Charger for WiCo, USB A-Micro B Cable, USB A-B Cable, Screwdriver
- NBR/PTFE로 구성된 Diaphragm으로 내화학성 우수
- Control Type : 무선 컨트롤러(WiCo)를 통해 작동 가능
- VSE1(Optional)을 사용하여 Solvent가 실내 공기로 방출되지 않게 보호
- USB Interface, RS232 Interface
- Protection class according to DIN EN 60529 (IP20) • IKA

Type	Digital	Control
Pumping Speed	22 L/min	
Ultimate Pressure	2 mbar	
Inlet/Outlet	8mm	
Diaphragm Material	NBR/PTFE	
Diaphragm Valve Material	FFPM	
Noise Level	54 dB(A)	
Size (w×d×h)/Weight	150×370×375mm/11.5 kg	270×500×435 mm/14 kg
Environmental Conditions	Temp. 5-40°C / Humidity 80%	
Power Input	130W	140W
Power	100-240V, 50/60Hz	

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
V08-05-031	VC 10 lite	Vacuum Controller, Only w/ VACSTAR Digital	1,900,000
P06-05-013	20015406	VSW 1 Vacuum Safety Wouff Bottle	910,000
P06-05-015	35001873	VSE 1 Vacuum Safety Emission Condenser	1,060,000
P06-05-025	35002702	VC 10.500 Stand pole(For mounting the VC 10 lite to the VACSTAR digital)	130,000



20005194



20004850



20004893



20013886



Digital type

Control type



VC 10 lite



20015406



35001873