

Vacuum Pump

표준형 진공 펌프, Edwards

Cat. No.	Model	Pumping Speed	Motor Power	EA Price ₩
P06-121-103	ED-RV3	65 L/min	550W	3,421,700
P06-121-105	ED-RV5	103 L/min		3,506,000
P06-121-108	ED-RV8	167 L/min		3,783,100
P06-121-112	ED-RV12	238 L/min		4,048,200

- 내화학성 Full Teflon 재질로, 부식성이 강한 유체 사용에 이상적임
- 용도 : 1 Lit. 정도의 소량 Filtration용
- Water Jet Pump나 Aspirator 대체용
- Rotary Evaporator와 같이 사용하기에 이상적임
- 환경 친화적이며, 증기와 응축에 강함
- 2 Stage (단, N810FT, 820FT, 840FT 1 Stage), 220V, 60Hz
- Edwards



Cat. No.	P06-121-103	P06-121-105	P06-121-108	P06-121-112
Model	ED-RV3	ED-RV5	ED-RV8	ED-RV12
Pumping Speed	65 L/min	103 L/min	167 L/min	238 L/min
Ultimate Pressure	1.5×0.001 Torr	1.5×0.001 Torr	1.5×0.001 Torr	1.5×0.001 Torr
Ultimate Pressure GB1	2.3×0.01 Torr	2.3×0.01 Torr	2.3×0.01 Torr	2.3×0.01 Torr
Ultimate Pressure GB2	9.1×0.01 Torr	9.1 X 0.01 Torr	4.6×0.01 Torr	9.1 X 0.01 Torr
Motor Power	550 W			
Motor Speed	1800 rpm			
Oil Capacity	700 ml		750 ml	1 L
Inlet	NW25, Nozzle 12 mm (기본 포함)			
Outlet	NW25			
무게	25 kg		28 kg	29 kg
Operation Temp. Range	12 ~ 40 °C			
Size	w158×d430×h221 mm		w158×d470×h221 mm	w158×d490×h221 mm
Voltage	220V, 60Hz			

Vacuum Pump

대용량 진공 펌프, Edwards

Cat. No.	Model	Pumping Speed	Motor Power	EA Price ₩
P06-121-201	ED-E2M18	340 L/min	750W	5,361,400
P06-121-211	ED-E2M28	550 L/min		5,686,700

- 대용량 진공 펌프
- 역류 방지 장치 (Check Valve)가 부착되어 있음
- Gas Ballast Valve를 가지고 있음
- 2 Stage Rotary Vane Pump
- Power Cable 포함
- 공통 약세사리는 "Edwards Accessory"란참조
- Edwards



Cat. No.	P06-121-201	P06-121-211
Model	ED-E2M18	ED-E2M28
Pumping Speed	340 L/min	550 L/min
Ultimate Pressure - without Gas Ballast - with Gas Ballast	7.4×0.0001 Torr 1.1×0.01 Torr	7.5×0.0001 Torr 1.1×0.01 Torr
Motor Power	750 W	
Motor Speed	1720 rpm	
Oil Capacity	1050 ml	1500 ml
Inlet	NW25, Nozzle 15 mm (기본 장착)	
Outlet	3/8" BSP Tapped Hole	
무게	39 kg	40 kg
Operation Temp. Range	12 ~ 40°C	
Size (w×d×h)	171×550×251 mm	160×584×241 mm
Voltage	220V, 60Hz	