

Vacuum Pump 진공 펌프, Vacuumer®

Cat. No.	Model	Pumping Speed	Motor Power	EA Price ₩
P06-62-702	VOP-50	50 L/min	200 W	1,450,000
P06-62-706	VOP-60	60 L/min		1,550,000
P06-62-710	VOP-100	120 L/min	400 W	1,650,000
P06-62-715	VOP-150	200 L/min	550 W	1,900,000
P06-62-720	VOP-200	260 L/min		2,000,000

- 고급형 고진공
- Gas Ballast Valve 부착
- Double Geage Vacuum Pump
- SUS 손잡이 장착
- Vacuumer®
- 과부하 보호장치, 역류방지장치 내장
- 저소음 저진동형, 강제공냉, 급유방식
- Thermal Protector Equiped
- 슬러지 오염방지 필터 장착



VOP-50



VOP-150

Oil Mist Trap for VOP, NW25

Cat. No.	Model	적용 기종	규격	EA Price ₩
P06-62-815	OMT-150P	VOP-60, VOP-100	200 L/min 이하	289,200
P06-62-820	OMT-300P	VOP-150, VOP-200	560 L/min 이하	433,700

- NW25 규격의 Oil Mist Trap 사용 가능



P06-62-815

P06-62-820

Cat. No.	P06-62-702	P06-62-706	P06-62-710	P06-62-715	P06-62-720
Model	VOP-50	VOP-60	VOP-100	VOP-150	VOP-200
Pumping Speed(50/60Hz)	42/50 L/min	50/60 L/min	100/120 L/min	170/200 L/min	220/260 L/min
Ultimate Pressure	6.7×10^{-2} (6.7×10^{-1}) Pa				
Motor Power	200 W		400 W	550 W	
Motor Speed (50/60 Hz)	1430 / 1730 rpm				
Oil Capacity	800 ml		900 ml	1 L	1.1 L
Inlet / Outlet	NW25 / Hose Nipple od12 mm 포함			NW25 / Hose Nipple od15 mm 포함	
Temp. Range	7 ~ 40 °C				
Weight	17 kg	17 kg	25 kg	28.5 kg	30 kg
Size (w × d × h, mm)	170 × 398 × 235	165 × 411 × 234	171 × 464 × 238	189 × 504 × 253	189 × 524 × 253
Voltage	220V, 50/60Hz				

Vacuum Pump 보급형 진공 펌프

Cat. No.	Model	Pumping Speed	Power	EA Price ₩
P06-197-100	W2V10	100 L/min	0.4 kW	1,132,500
P06-197-110	W2V20	200 L/min		1,289,200

- 국산 경제형 진공 펌프
- 구조가 간단하며 호환성이 양호
- 확실한 진공 설계와 잔류 공기 배출 (Gas Ballast) 설계
- OD12mm Hose Nipple과 Clamp 포함
- Oil Mist Trap은 별도 구매
- 400L/min ~ 600L/min 제품은 홈페이지 참조
- LK Labkorea
- Compact 사이즈 및 소음, 진동이 작음



EMF10

Cat. No.	P06-197-100	P06-197-110
Model	W2V10	W2V20
Pumping Speed	100 L/min	200 L/min
진공도	$\leq 1 \times 10^{-3}$ Torr (6.7×10^{-1} Pa), Gas Ballast Closed / 5×10^{-2} Torr (6.7 Pa), Gas Ballast Opened	
Full Load Power	0.4 kW (0.5 HP)	
Motor Speed	1700 rpm	
Oil Capacity	500 ml	600 ml
Weight	19 kg	21 kg
Inlet / Outlet	NW25 / NW25	
사용 온도 조건	7 ~ 40 °C	
Size (mm)	w150 x d398 x h251	w150 x d434 x h251
Power	220V, 50/60Hz, 1Ø	



WOF150

Optional Oil Mist Trap

Cat. No.	Model	규격	Size	EA Price ₩
P06-197-612	WOF150	NW25 용	Ø90 x h148 mm	120,500
P06-198-512	EMF10	200L 이하	Ø97 x h171 mm	289,200