

## Wire Washing Bottle Rack / Carrier 스텐션 워싱 바틀 랙

Cat. No.	Model	Hole	Size (w×d×h)	EA Price ₩
R01-11-112	WBR-S6	Ø75 mm, 6홀	290×235×185 mm	34,400
R01-11-114	WBR-S4	Ø95 mm, 4홀	220×220×110 mm	35,000

- 재질 : Stainless Steel
- WBR-S6은 500 ml 세척병의 운반, 보관에 알맞음
- WBR-S4은 1 L 세척병의 운반, 보관에 알맞음
- LK Labkorea



## Wire Rack for B.O.D Bottle BOD 바틀용 와이어 랙

Cat. No.	Model	사용 BOD Bottle	Size (w×d×h)	EA Price ₩
B08-53-111	LK06020	60 ml × 20 개	340×220×150 mm	55,700
B08-55-114	LK30012	300 ml × 12 개	340×250×195 mm	40,200

- 재질 : Stainless Steel
- LK06020 : up to Ø45mm
- LK30012 : up to Ø71mm
- LK Labkorea



## Utility Carrier / Bottle Rack 다용도 운반대

Cat. No.	Model	내부 Size	EA Price ₩
R01-46-111	T240080	w354×d100×h110 mm (2열)	36,900

- 재질 : PP, -10 ~ 135 °C 내열, Autoclavable
- 각종 기구 및 병 등의 운반 및 보관용
- 전체 Size : w380×d240×h115 mm
- Tarsons



## Safety Solvent Bottle Carrier 안전 솔벤트 바틀 캐리어, Tote Style

Cat. No.	Model	Type	색상	EA Price ₩
B08-134-970	235349	500 ml & 1 L 용	Red	135,200
B08-134-974	235350		Black	135,200
B08-134-978	235345	2.5 L & 4 L 용	Red	216,900
B08-134-982	235346		White	216,900
B08-134-986	235348		Black	216,900

- 재질 : Thermoplastic Elastomer, -45 ~ 110°C 내열, Autoclavable
- 자세한 특징은 169 page 참조
- Qorpak



## Safety Solvent Bottle Carrier 안전 솔벤트 바틀 캐리어, Cover Style

- 재질 : Epoxy Coated PE, Heavy-gauge Steel Wire Bail
- Wide Ribs in Both the Sides and on The Bottom to Support and Cushion Bottles
- Cone Style Lid : 16mm Opening in the Top for Insertion of Tubing - Simplifies the Transfer of Liquids without Pouring
- Truncated Style Lid : Pouring without the Removal of the Lid
- Bel-Art

### Cone Style, 16 mm Opening for Tubing Inserting ①

Cat. No.	Model	적용 Bottle	Size	EA Price ₩
B08-11-430	F16957-0000	2.5 L	Ø180×h350 mm	260,900
B08-11-434	F16958-0000	4 L	Ø200×h370 mm	309,500

### Truncated Style, for Pouring ②

Cat. No.	Model	적용 Bottle	Size	EA Price ₩
B08-11-438	F16959-0000	2.5 L	Ø180×h320 mm	260,900
B08-11-442	F16960-0000	4 L	Ø200×h350 mm	307,800

