

Glove Holder / Glove Dispenser 글러브 홀더 / 글러브 디스펜서

Cat. No.	Model	Holder	Size (mm)	EA Price ₩
R01-11-178	370050	1 Box	w260×d105×h140	37,800
R01-11-181	370060	2 Box	w260×d105×h265	52,000
R01-11-184	370070	3 Box	w260×d105×h400	66,100

- 재질 : 투명 Acrylic
- 벽에 걸어 사용할 수 있음
- LK Labkorea



Separatory Funnel Rack 분액 여두 원형 랙 ①

Cat. No.	Model	분액 여두	홀 규격	EA Price ₩
R01-143-088	93-3003-5	50 ml 용	Ø45 mm × 6개	99,900
R01-143-090	93-3006-5		Ø45 mm × 10개	135,300
R01-143-096	93-3004-5	100 / 250 ml 용	Ø60 mm × 6개	129,200
R01-143-098	93-3007-5		Ø60 mm × 10개	162,400
R01-143-102	93-3005-5	500 ml / 1 L 용	Ø83 mm × 6개	166,500

- 재질 : PVC 하부 다이, MF 상부 거치대
- 안정된 Base와 높이 조절 가능한 지지판, 견고한 기둥으로 구성
- Base Size : Ø300 mm
(단, 93-3007-5 : Ø340 mm, 93-3005-5 : Ø320 mm)
- Height : 530 mm
- Sansyo



Funnel Stand 여두용 스탠드 ②

Cat. No.	Model	홀 Size	Size (w×d×h)	EA Price ₩
R01-143-303	93-3001-5	Ø44 mm	240×90×428 mm	39,500

- 재질 : PVC
- Sansyo

Separatory Funnel Rack / Funnel Rack 스텐션 여두 랙 / 여두대

Cat. No.	Model	홀 Size	홀 수	높이	EA Price ₩
R01-11-190	SFR-06S	Ø70 mm	6개 ①	190 ~ 300 mm	55,000
R01-11-196	SFR-06M	Ø90 mm		190 ~ 300 mm	55,000
R01-11-198	SFR-06L	Ø115 mm		300 ~ 500 mm	145,800
R01-11-200	SFR-06XL	Ø140 mm		300 ~ 500 mm	173,300
R01-11-202	SFR-12S	Ø70 mm	12개 ②	190 ~ 300 mm	77,500
R01-11-208	SFR-12M	Ø90 mm		190 ~ 300 mm	77,500
R01-11-210	SFR-12L	Ø115 mm		300 ~ 500 mm	207,400
R01-11-212	SFR-12XL	Ø140 mm		300 ~ 500 mm	216,700

- 재질 : Stainless Steel
- 양 옆의 나사를 이용하여 높이를 조절함
- LK Labkorea



Parafilm Dispenser 파라필름 디스펜서

Cat. No.	Model	재질	Size (w×d×h)	EA Price ₩
W06-11-043	PD-171	ABS	125×155×146 mm	32,500
W06-11-047	PD-171B	Replacement Blade		5,100

- 재질 : ABS (Acrylonitrile Butadiene Styrene Copolymer)
- 칼날이 안쪽에 숨어있어 파라필름 절단할 때 안전함
- LK Labkorea

