

## Test Tube Rack with Handle 핸들형 테스트 튜브 랙

Cat. No.	Model	홀 (Type)	색상	EA Price ₩	PK Price ₩
R01-46-401	T202010-B	31홀	Blue	8,200	4ea 31,100
R01-46-403	T202010-W	Ø13mm용	White	8,200	4ea 31,100
R01-46-411	T202020-B	31홀	Blue	8,200	4ea 31,100
R01-46-413	T202020-W	Ø16mm용	White	8,200	4ea 31,100
R01-46-421	T202030-B	20홀	Blue	10,800	4ea 41,200
R01-46-423	T202030-W	Ø20mm용	White	10,800	4ea 41,200
R01-46-431	T202040-B	18홀	Blue	10,800	4ea 41,200
R01-46-433	T202040-W	Ø25mm용	White	10,800	4ea 41,200
R01-46-441	T202050-B	12홀	Blue	10,800	4ea 41,200
R01-46-444	T202050-W	Ø32mm용	White	10,800	4ea 41,200



- 재질 : PP, Autoclavable
- 손잡이가 있어 이동이 편리함
- Size : (31홀) w272×d65×h110 mm  
(20홀, 18홀) w322×d69×h129 mm  
(12홀) w278×d83×h126 mm
- Tarsons



## Plastic Test Tube Rack 플라스틱 튜브 랙


Cat. No.	Model	홀 (Type)	색상	EA Price ₩	PK Price ₩
R01-20-707	4340001	84홀 (A) (6×14홀) Ø13 mm용	Blue	21,200	5ea 98,800
R01-20-713	4340000		White	21,200	5ea 98,800
R01-20-715	4340002		Red	21,200	5ea 98,800
R01-20-781	4340061	55홀 (A) (5×11홀) Ø16 mm용	Blue	21,200	5ea 98,800
R01-20-783	4340060		White	21,200	5ea 98,800
R01-20-785	4340062		Red	21,200	5ea 98,800
R01-20-719	4340011	55홀 (A) (5×11홀) Ø18 mm용	Blue	21,200	5ea 98,800
R01-20-725	4340010		White	21,200	5ea 98,800
R01-20-727	4340012		Red	21,200	5ea 98,800
R01-20-731	4340021	40홀 (A) (4×10홀) Ø20 mm용	Blue	21,200	5ea 98,800
R01-20-747	4340020		White	21,200	5ea 98,800
R01-20-749	4340022		Red	21,200	5ea 98,800
R01-20-753	4340031	32홀 (B) (4×8홀) Ø25 mm용	Blue	21,200	5ea 98,800
R01-20-759	4340030		White	21,200	5ea 98,800
R01-20-761	4340032		Red	21,200	5ea 98,800
R01-20-765	4340041	21홀 (B) (3×7홀) Ø30 mm용	Blue	21,200	5ea 98,800
R01-20-771	4340040		White	21,200	5ea 98,800
R01-20-773	4340042		Red	21,200	5ea 98,800

- 재질 : PP, -10 ~ 135 °C 내열, Autoclavable
- 쌓아 올릴 수 있으며, 문자, 숫자 식 위치 표시
- 물에 넣어도 뜨지 않음
- Size : w265×d126 mm, Height : 75 mm (Type A), 88 mm (Type B)
- Brand GmBH



## Polygrid® Test Tube Rack 테스트 튜브 랙

Cat. No.	Model	홀 (Type)	색상	EA Price ₩	PK Price ₩
R01-46-030	T205010-W	90홀 Ø13mm용	White	11,300	8ea 86,300
R01-46-033	T205010-B		Blue	11,300	8ea 86,300
R01-46-036	T205010-R		Red	11,300	8ea 86,300
R01-46-039	T205010-Y		Yellow	11,300	8ea 86,300
R01-46-042	T205020-W	60홀 Ø16mm용	White	11,300	8ea 86,300
R01-46-045	T205020-B		Blue	11,300	8ea 86,300
R01-46-048	T205020-R		Red	11,300	8ea 86,300
R01-46-051	T205020-Y		Yellow	11,300	8ea 86,300
R01-46-054	T205030-W	40홀 Ø20mm용	White	11,300	8ea 86,300
R01-46-057	T205030-B		Blue	11,300	8ea 86,300
R01-46-060	T205030-R		Red	11,300	8ea 86,300
R01-46-063	T205030-Y		Yellow	11,300	8ea 86,300
R01-46-066	T205040-W	24홀 Ø25mm용	White	11,300	8ea 86,300
R01-46-069	T205040-B		Blue	11,300	8ea 86,300
R01-46-072	T205040-R		Red	11,300	8ea 86,300
R01-46-075	T205040-Y		Yellow	11,300	8ea 86,300



- 재질 : RPP(Re-Enforced PP), 0 ~ 135 °C 내열 Autoclavable
- Size : w246×d104×h64 mm
- Tarsons