

I Reaction Vessel

비커형 반응조, LukeGL[®]

- Borosilicate Glass 3.3
- 비커 Type의 환저 또는 평저형 반응조
- O-ring 및 Clamp는 별도 구매해야 함 (Accessory란 참조)
- Cover의 Neck Size는 선택 가능함
- 500 ml ~ 5 L 는 용량 표시용 눈금을 가지고 있음
- 상부 Flange에 O-ring홈이 있는것과 없는것 2 종류가 있음
- LukeGL[®]

일반 비커형 평저 반응조 ①

Cat. No.	Model	용량	외경	전체높이	Flange	EA Price ₩
R03-09-007	GRG1003	300 ml	85 mm	90 mm	115 mm	48,900
R03-09-014	GRG1005	500 ml	90 mm	130 mm	120 mm	66,300
R03-09-021	GRG1010	1 L	110 mm	150 mm	145 mm	79,500
R03-09-028	GRG1020	2 L	130 mm	200 mm	170 mm	121,400
R03-09-035	GRG1030	3 L	150 mm	220 mm	190 mm	163,600
R03-09-042	GRG1050	5 L	180 mm	260 mm	220 mm	258,100
R03-09-049	GRG1070	7 L	180 mm	350 mm	220 mm	337,300
R03-09-056	GRG1100	10 L	210 mm	360 mm	250 mm	626,500
R03-09-063	GRG1200	20 L	250 mm	530 mm	290 mm	1,012,000

오링 비커형 평저 반응조 ②

Cat. No.	Model	용량	외경	전체높이	Flange	EA Price ₩
R03-09-070	GRO1003	300 ml	85 mm	95 mm	115 mm	75,900
R03-09-077	GRO1005	500 ml	90 mm	135 mm	120 mm	101,200
R03-09-084	GRO1010	1 L	110 mm	155 mm	140 mm	114,500
R03-09-091	GRO1020	2 L	130 mm	205 mm	160 mm	156,600
R03-09-098	GRO1030	3 L	150 mm	225 mm	180 mm	209,600
R03-09-105	GRO1050	5 L	180 mm	265 mm	210 mm	307,200
R03-09-112	GRO1070	7 L	180 mm	355 mm	210 mm	481,900
R03-09-119	GRO1100	10 L	210 mm	365 mm	240 mm	632,500
R03-09-126	GRO1200	20 L	250 mm	535 mm	290 mm	1,097,600

일반 비커형 환저 반응조 ③

Cat. No.	Model	용량	외경	전체높이	Flange	EA Price ₩
R03-09-133	GRG2003	300 ml	85 mm	110 mm	115 mm	49,400
R03-09-140	GRG2005	500 ml	90 mm	150 mm	120 mm	77,100
R03-09-147	GRG2010	1 L	110 mm	170 mm	145 mm	92,800
R03-09-154	GRG2020	2 L	130 mm	220 mm	170 mm	121,700
R03-09-161	GRG2030	3 L	150 mm	250 mm	190 mm	163,600
R03-09-168	GRG2050	5 L	180 mm	300 mm	220 mm	291,600
R03-09-175	GRG2070	7 L	180 mm	380 mm	220 mm	395,200
R03-09-182	GRG2100	10 L	210 mm	375 mm	250 mm	715,700
R03-09-189	GRG2200	20 L	250 mm	560 mm	290 mm	1,041,000

오링 비커형 환저 반응조 ④

Cat. No.	Model	용량	외경	전체높이	Flange	EA Price ₩
R03-09-196	GRO2003	300 ml	85 mm	115 mm	115 mm	75,900
R03-09-203	GRO2005	500 ml	90 mm	155 mm	120 mm	101,200
R03-09-210	GRO2010	1 L	110 mm	175 mm	140 mm	114,500
R03-09-217	GRO2020	2 L	130 mm	225 mm	160 mm	156,600
R03-09-224	GRO2030	3 L	150 mm	255 mm	180 mm	209,600
R03-09-231	GRO2050	5 L	180 mm	305 mm	210 mm	307,200
R03-09-238	GRO2070	7 L	180 mm	385 mm	210 mm	372,300
R03-09-245	GRO2100	10 L	210 mm	380 mm	240 mm	632,500
R03-09-252	GRO2200	20 L	250 mm	565 mm	290 mm	1,097,600

