

## Baffled Reaction Vessel 비커형 베플 반응조, LukeGL<sup>®</sup>

- Borosilicate Glass 3.3
- 비커 Type의 Baffle형 환저 또는 평저 반응조
- O-ring 및 Clamp는 별도 구매해야 함 (Accessory란 참조)
- Cover의 Neck Size는 선택 가능함
- 상부 Flange에 O-ring홈이 있는것과 없는것 2종류가 있음
- LukeGL<sup>®</sup>

### 일반 비커형 평저 베플 반응조 ①

Cat. No.	Model	용량	외경	전체높이	Flange	EA Price ₩
R03-09-259	GRG3003	300 ml	85 mm	90 mm	115 mm	71,100
R03-09-266	GRG3005	500 ml	90 mm	130 mm	120 mm	116,400
R03-09-273	GRG3010	1 L	110 mm	150 mm	145 mm	128,200
R03-09-280	GRG3020	2 L	130 mm	200 mm	170 mm	165,300
R03-09-287	GRG3030	3 L	150 mm	220 mm	190 mm	203,600
R03-09-294	GRG3050	5 L	180 mm	260 mm	220 mm	313,300
R03-09-301	GRG3070	7 L	180 mm	350 mm	220 mm	481,900
R03-09-308	GRG3100	10 L	210 mm	360 mm	250 mm	715,700
R03-09-315	GRG3200	20 L	250 mm	530 mm	290 mm	1,151,800

### 오링 비커형 평저 베플 반응조 ②

Cat. No.	Model	용량	외경	전체높이	Flange	EA Price ₩
R03-09-322	GRO3003	300 ml	85 mm	95 mm	115 mm	101,800
R03-09-329	GRO3005	500 ml	90 mm	135 mm	120 mm	144,100
R03-09-336	GRO3010	1 L	110 mm	155 mm	140 mm	156,600
R03-09-343	GRO3020	2 L	130 mm	205 mm	160 mm	188,000
R03-09-350	GRO3030	3 L	150 mm	225 mm	180 mm	234,900
R03-09-357	GRO3050	5 L	180 mm	265 mm	210 mm	379,000
R03-09-364	GRO3070	7 L	180 mm	355 mm	210 mm	524,700
R03-09-371	GRO3100	10 L	210 mm	365 mm	240 mm	772,300
R03-09-378	GRO3200	20 L	250 mm	535 mm	290 mm	1,159,000

### 일반 비커형 환저 베플 반응조 ③

Cat. No.	Model	용량	외경	전체높이	Flange	EA Price ₩
R03-09-385	GRG4003	300 ml	85 mm	110 mm	115 mm	71,100
R03-09-392	GRG4005	500 ml	90 mm	150 mm	120 mm	121,700
R03-09-399	GRG4010	1 L	110 mm	170 mm	145 mm	134,900
R03-09-406	GRG4020	2 L	130 mm	220 mm	170 mm	173,500
R03-09-413	GRG4030	3 L	150 mm	250 mm	190 mm	203,600
R03-09-420	GRG4050	5 L	180 mm	300 mm	220 mm	327,700
R03-09-427	GRG4070	7 L	180 mm	380 mm	220 mm	481,900
R03-09-434	GRG4100	10 L	210 mm	380 mm	250 mm	726,500
R03-09-441	GRG4200	20 L	250 mm	560 mm	290 mm	1,151,800

### 오링 비커형 환저 베플 반응조 ④

Cat. No.	Model	용량	외경	전체높이	Flange	EA Price ₩
R03-09-448	GRO4003	300 ml	85 mm	115 mm	115 mm	101,800
R03-09-455	GRO4005	500 ml	90 mm	155 mm	120 mm	144,100
R03-09-462	GRO4010	1 L	110 mm	175 mm	140 mm	156,600
R03-09-469	GRO4020	2 L	130 mm	225 mm	160 mm	188,000
R03-09-476	GRO4030	3 L	150 mm	255 mm	180 mm	234,900
R03-09-483	GRO4050	5 L	180 mm	305 mm	210 mm	379,000
R03-09-490	GRO4070	7 L	180 mm	385 mm	210 mm	524,700
R03-09-497	GRO4100	10 L	210 mm	385 mm	240 mm	772,300
R03-09-504	GRO4200	20 L	250 mm	565 mm	290 mm	1,127,700

