

## Reaction Vessel, Bottle type

### 바틀형 반응조, 10 L & 20 L, LukeGL<sup>®</sup>


- Borosilicate Glass 3.3
- O-ring 및 Clamp는 별도임 (Accessory란 참조), Cover Neck Size 선택가능
- LukeGL<sup>®</sup>

#### 일반형 여과병 반응조

Cat. No.	Model	용량	외경	전체높이	Flange	EA Price ₩
R03-09-507	GRG5100	10 L	240 mm	360 mm	190 mm	619,600
R03-09-511	GRG5200	20 L	299 mm	450 mm	220 mm	1,015,500

#### 오링형 여과병 반응조

Cat. No.	Model	용량	외경	전체높이	Flange	EA Price ₩
R03-09-515	GRO5100	10 L	240 mm	365 mm	180 mm	671,200
R03-09-520	GRO5200	20 L	299 mm	455 mm	210 mm	1,084,300



오링형

## Reaction Vessel, Flask type

### 플라스크형 반응조, LukeGL<sup>®</sup>

- Borosilicate Glass 3.3
- O-ring 및 Clamp는 별도임 (Accessory란 참조), Cover Neck Size 선택가능
- LukeGL<sup>®</sup>

#### 일반형 플라스크 반응조 ①

Cat. No.	Model	용량	외경	전체높이	Flange	EA Price ₩
R03-09-525	GRG6005	500 ml	105 mm	130 mm	115 mm	58,000
R03-09-532	GRG6010	1 L	131 mm	155 mm	135 mm	92,800
R03-09-539	GRG6020	2 L	166 mm	185 mm	145 mm	125,300
R03-09-546	GRG6030	3 L	185 mm	210 mm	160 mm	178,300
R03-09-553	GRG6070	5 L	223 mm	240 mm	170 mm	242,500
R03-09-560	GRG6100	10 L	279 mm	310 mm	190 mm	524,100
R03-09-567	GRG6200	20 L	345 mm	370 mm	220 mm	744,000

#### 오링형 플라스크 반응조 ②

Cat. No.	Model	용량	외경	전체높이	Flange	EA Price ₩
R03-09-574	GRO6005	500 ml	105 mm	130 mm	115 mm	92,800
R03-09-581	GRO6010	1 L	131 mm	155 mm	130 mm	113,300
R03-09-588	GRO6020	2 L	166 mm	185 mm	140 mm	156,600
R03-09-595	GRO6030	3 L	185 mm	210 mm	150 mm	225,300
R03-09-602	GRO6070	5 L	223 mm	240 mm	160 mm	286,100
R03-09-609	GRO6100	10 L	279 mm	310 mm	180 mm	572,300
R03-09-616	GRO6200	20 L	345 mm	370 mm	210 mm	795,200



## Reaction Vessel with Needle Valve, Flask type

### 오링형 밸브식 플라스크형 반응조, LukeGL<sup>®</sup>

Cat. No.	Model	용량	외경	전체높이	Flange	EA Price ₩
R03-09-924	GRO9005	500 ml	109 mm	130 mm	115 mm	235,700
R03-09-931	GRO9010	1 L	137 mm	155 mm	130 mm	256,000
R03-09-938	GRO9020	2 L	161 mm	185 mm	140 mm	295,200
R03-09-945	GRO9030	3 L	184 mm	210 mm	150 mm	348,000
R03-09-952	GRO9050	5 L	223 mm	240 mm	160 mm	472,900
R03-09-959	GRO9100	10 L	279 mm	310 mm	180 mm	843,400
R03-09-966	GRO9200	20 L	345 mm	370 mm	210 mm	1,067,800

- Borosilicate Glass 3.3
- O-ring 및 Clamp는 별도 구매해야 함 (Accessory란 참조)
- 높이는 밸브의 높이를 제외한 전체 높이임
- 요청에 의해 일반형 Flange 제품도 공급가능함, Cover Neck Size 선택가능
- LukeGL<sup>®</sup>

