

## Jacketed Reaction Vessel 자켓식 반응조, LukeGL<sup>®</sup>

- Borosilicate Glass 3.3
- O-ring 및 Clamp는 별도 구매해야 함 (Accessory란 참조)
- 500 ml ~ 5 L는 용량 표시용 눈금을 가지고 있음
- Connection : GL14 (100 ml ~ 1 L), GL18 (2 L ~ 20 L)
- LukeGL<sup>®</sup>

### 일반형 이중 자켓식 평저 반응조 ①

Cat. No.	Model	용량	외경 (mm)		높이 (mm)		Flange (mm)	EA Price ₩
			겉관	속관	겉관	속관		
R04-09-006	JRG1001	100 ml	85	60	110	85	90	177,100
R04-09-012	JRG1002	250 ml	110	80	140	105	110	177,100
R04-09-018	JRG1005	500 ml	120	90	160	110	120	199,900
R04-09-024	JRG1010	1 L	140	110	200	140	145	244,600
R04-09-030	JRG1020	2 L	160	130	250	190	170	359,300
R04-09-036	JRG1030	3 L	180	150	270	210	190	492,500
R04-09-042	JRG1050	5 L	200	170	340	270	210	728,900
R04-09-048	JRG1070	7 L	200	170	440	370	210	974,700
R04-09-054	JRG1100	10 L	250	210	400	330	250	1,351,800
R04-09-060	JRG1200	20 L	300	250	560	480	290	3,108,400

### 오링형 이중 자켓식 평저 반응조 ②

Cat. No.	Model	용량	외경 (mm)		높이 (mm)		Flange (mm)	EA Price ₩
			겉관	속관	겉관	속관		
R04-09-066	JRO1001	100 ml	85	60	110	70	90	177,100
R04-09-072	JRO1002	250 ml	110	80	140	85	110	203,600
R04-09-078	JRO1005	500 ml	120	90	165	110	120	244,600
R04-09-084	JRO1010	1 L	140	110	205	140	140	269,900
R04-09-090	JRO1020	2 L	160	130	255	190	160	409,600
R04-09-096	JRO1030	3 L	180	150	275	210	180	561,400
R04-09-102	JRO1050	5 L	200	170	345	270	200	814,500
R04-09-108	JRO1070	7 L	200	170	445	370	200	939,800
R04-09-114	JRO1100	10 L	250	210	405	330	240	1,506,000
R04-09-120	JRO1200	20 L	300	250	565	480	290	3,100,000

### 일반형 이중 자켓식 환저 반응조 ③

Cat. No.	Model	용량	외경 (mm)		높이 (mm)		Flange (mm)	EA Price ₩
			겉관	속관	겉관	속관		
R04-09-126	JRG2001	100 ml	85	60	120	95	90	177,100
R04-09-132	JRG2002	250 ml	110	80	150	115	110	177,100
R04-09-138	JRG2005	500 ml	120	90	175	125	120	209,600
R04-09-144	JRG2010	1 L	140	110	210	155	145	256,600
R04-09-150	JRG2020	2 L	160	130	270	210	170	378,300
R04-09-156	JRG2030	3 L	180	150	285	230	190	518,100
R04-09-162	JRG2050	5 L	200	170	360	300	210	762,700
R04-09-168	JRG2070	7 L	200	170	475	400	210	860,200
R04-09-174	JRG2100	10 L	250	210	430	360	250	1,416,900
R04-09-180	JRG2200	20 L	300	250	630	520	290	3,010,800

### 오링형 이중 자켓식 환저 반응조 ④

Cat. No.	Model	용량	외경 (mm)		높이 (mm)		Flange (mm)	EA Price ₩
			겉관	속관	겉관	속관		
R04-09-186	JRO2001	100 ml	85	60	120	95	90	177,100
R04-09-192	JRO2002	250 ml	110	80	155	100	110	203,600
R04-09-198	JRO2005	500 ml	120	90	180	125	120	244,600
R04-09-204	JRO2010	1 L	140	110	215	155	140	269,900
R04-09-210	JRO2020	2 L	160	130	275	210	160	409,600
R04-09-216	JRO2030	3 L	180	150	290	230	180	561,400
R04-09-222	JRO2050	5 L	200	170	365	300	200	814,500
R04-09-228	JRO2070	7 L	200	170	480	400	200	956,600
R04-09-234	JRO2100	10 L	250	210	425	360	240	1,506,000
R04-09-240	JRO2200	20 L	300	250	630	520	290	3,100,000

