

Cover for Reaction Vessel

일반형 반응조 뚜껑, LukeGL®

- Borosilicate Glass 3.3
- 일반형 반응기 하부용 Cover
- Grease나 Gasket (Silicon or 테프론 재질)과 같이 사용
- 기본 Joint Size는 24/40임
- 기본 규격외의 Neck 종류, 모양, 규격은 요청에 의해 공급 가능
- LukeGL®



Flange		1-neck	2-neck	3-neck	4-neck	5-neck	6-neck
90 mm	Cat. No.	R04-09-H004	R04-09-H008	R04-09-H012	-	-	-
	Model	GRCG01085	GRCG02085	GRCG03085	-	-	-
	EA Price ₩	71,100	77,100	81,000			
110 mm	Cat. No.	R04-09-H028	R04-09-H032	R04-09-H036	R04-09-H040	-	-
	Model	GRCG01100	GRCG02100	GRCG03100	GRCG04100	-	-
	EA Price ₩	71,100	77,100	84,300	92,800		
115 mm	Cat. No.	R04-09-H052	R04-09-H056	R04-09-H060	R04-09-H064	-	-
	Model	GRCG01115	GRCG02115	GRCG03115	GRCG04115	-	-
	EA Price ₩	63,900	72,300	77,100	100,000		
120 mm	Cat. No.	R04-09-H076	R04-09-H080	R04-09-H084	R04-09-H088	R04-09-H092	-
	Model	GRCG01120	GRCG02120	GRCG03120	GRCG04120	GRCG05120	-
	EA Price ₩	71,100	77,100	84,300	94,000	106,000	
135 mm	Cat. No.	R04-09-H100	R04-09-H104	R04-09-H108	R04-09-H112	R04-09-H116	R04-09-H120
	Model	GRCG01135	GRCG02135	GRCG03135	GRCG04135	GRCG05135	GRCG06135
	EA Price ₩	70,900	83,100	95,200	100,000	109,600	119,300
145 mm	Cat. No.	R04-09-H124	R04-09-H128	R04-09-H132	R04-09-H136	R04-09-H140	R04-09-H144
	Model	GRCG01145	GRCG02145	GRCG03145	GRCG04145	GRCG05145	GRCG06145
	EA Price ₩	73,600	85,500	97,600	104,800	113,300	126,500
160 mm	Cat. No.	R04-09-H148	R04-09-H152	R04-09-H156	R04-09-H160	R04-09-H164	R04-09-H168
	Model	GRCG01160	GRCG02160	GRCG03160	GRCG04160	GRCG05160	GRCG06160
	EA Price ₩	80,700	92,800	101,200	109,600	118,100	127,700
170 mm	Cat. No.	R04-09-H172	R04-09-H176	R04-09-H180	R04-09-H184	R04-09-H188	R04-09-H192
	Model	GRCG01170	GRCG02170	GRCG03170	GRCG04170	GRCG05170	GRCG06170
	EA Price ₩	89,400	97,900	108,400	118,100	130,100	138,600
190 mm	Cat. No.	R04-09-H196	R04-09-H200	R04-09-H204	R04-09-H208	R04-09-H212	R04-09-H216
	Model	GRCG01190	GRCG02190	GRCG03190	GRCG04190	GRCG05190	GRCG06190
	EA Price ₩	110,800	133,300	143,400	151,900	172,300	181,900
210 mm	Cat. No.	R04-09-H220	R04-09-H224	R04-09-H228	R04-09-H232	R04-09-H236	R04-09-H240
	Model	GRCG01210	GRCG02210	GRCG03210	GRCG04210	GRCG05210	GRCG06210
	EA Price ₩	160,200	169,900	180,700	190,400	202,400	212,000
220 mm	Cat. No.	R04-09-H244	R04-09-H248	R04-09-H252	R04-09-H256	R04-09-H260	R04-09-H264
	Model	GRCG01220	GRCG02220	GRCG03220	GRCG04220	GRCG05220	GRCG06220
	EA Price ₩	168,700	186,700	192,800	203,600	214,500	224,100
240 mm	Cat. No.	R04-09-H268	R04-09-H272	R04-09-H276	R04-09-H280	R04-09-H284	R04-09-H288
	Model	GRCG01240	GRCG02240	GRCG03240	GRCG04240	GRCG05240	GRCG06240
	EA Price ₩	268,700	278,300	288,000	298,800	310,800	320,500
250 mm	Cat. No.	R04-09-H292	R04-09-H296	R04-09-H300	R04-09-H304	R04-09-H308	R04-09-H312
	Model	GRCG01250	GRCG02250	GRCG03250	GRCG04250	GRCG05250	GRCG06250
	EA Price ₩	337,300	347,000	357,800	367,500	379,500	389,200
290 mm	Cat. No.	R04-09-H316	R04-09-H320	R04-09-H324	R04-09-H328	R04-09-H332	R04-09-H336
	Model	GRCG01290	GRCG02290	GRCG03290	GRCG04290	GRCG05290	GRCG06290
	EA Price ₩	675,900	686,700	698,800	708,400	718,100	730,100

조인트 규격

Joint Size							
10/30	14/20	24/40	19/38	29/32	29/42	34/45	45/50

Screw Connector Size

GL14	Connection : Ø6 ~ 8 mm Tubing용
GL18	Connection : Ø8 ~ 10 mm Tubing용

- 표준형과 고온용이 있음 (286 Page 참조)

