

O-ring for Reaction Vessel

반응조 오링



Viton O-ring

Cat. No.	Model	사용 Flange	EA Price ₩
R04-243-007	VO-F85	90	5,200
R04-243-014	VO-F100	110	9,800
R04-243-021	VO-F115	115	6,800
R04-243-028	VO-F120	120	10,000
R04-243-035	VO-F130	130	10,500
R04-243-042	VO-F140	140	10,800
R04-243-049	VO-F150	150	11,000
R04-243-056	VO-F160	160	11,500
R04-243-063	VO-F180	180	15,000
R04-243-070	VO-F200	200	23,200
R04-243-077	VO-F210	210	27,000
R04-243-084	VO-F230	230	29,000
R04-243-091	VO-F240	240	30,900
R04-243-098	115520	290	40,500

- Viton 재질 : 불소 고무 재질의 O-ring
- 사용온도 : -10℃ ~ 200℃
- 내열성, 내화학성이 뛰어난 불소 고무 오링
- 대부분의 화학약품 및 부식성 유체에 사용 가능함



Encapsulated Teflon O-ring

Cat. No.	Model	사용 Flange	EA Price ₩
R04-243-507	TO-F85	90	23,700
R04-243-514	TO-F110	110	28,500
R04-243-521	TO-F115	115	28,500
R04-243-528	TO-F120	120	30,900
R04-243-535	TO-F130	130	30,800
R04-243-542	TO-F140	140	33,200
R04-243-549	TO-F150	150	35,600
R04-243-556	TO-F160	160	38,200
R04-243-563	TO-F180	180	40,800
R04-243-570	TO-F200	200	46,000
R04-243-577	TO-F210	210	48,600
R04-243-584	TO-F230	230	51,400
R04-243-591	TO-F240	240	54,300
R04-243-598	TO-F290	290	62,700

- 테프론 Encapsulated O-ring
: Viton O-ring에 Teflon을 Capsule화 시킨 오링
- 사용온도 : -60℃ ~ 200℃

O-ring for Reaction Vessel

반응조 클램프

- Stainless Steel
- 간단하게 반응조 상, 하부를 고정 시킬수 있는 클램프임
- 간격 조절용 Screw Nut 포함



Cat. No.	Model	사용 Flange	EA Price ₩	Cat. No.	Model	사용 Flange	EA Price ₩
R04-09-912	SC-F90	90	30,900	R04-09-952	SC-F170	170	36,700
R04-09-916	SC-F110	110	32,800	R04-09-956	SC-F180	180	38,600
R04-09-920	SC-F115	115	32,800	R04-09-960	SC-F190	190	38,600
R04-09-924	SC-F120	120	32,800	R04-09-964	SC-F200	200	40,500
R04-09-928	SC-F130	130	34,700	R04-09-968	SC-F210	210	44,400
R04-09-932	SC-F135	135	34,700	R04-09-972	SC-F220	220	46,300
R04-09-936	SC-F140	140	34,700	R04-09-976	SC-F230	230	63,300
R04-09-940	SC-F145	145	36,700	R04-09-980	SC-F240	240	61,700
R04-09-944	SC-F150	150	36,700	R04-09-984	SC-F250	250	67,500
R04-09-948	SC-F160	160	36,700	R04-09-992	SC-F290	290	144,600