

Micro Reaction Vessel PFA 테프론 소형 반응기

Cat. No.	Model	용량	외부 Size	I.D.	EA Price ₩
R03-460-071	B318-40	90 ml	∅ 50 x h 96mm	∅ 45.5mm	1,084,700
R03-460-075	B318-64	240 ml	∅ 66 x h 130mm	∅ 60.0mm	1,276,100
R03-460-079	B318-80	500 ml	∅ 85 x h 158mm	∅ 80.0mm	1,990,400

- 재질 : PFA Body, PTFE Cover , -200 °C to 250°C 내열
- 테프론 재질로 최상의 내화학성을 가짐
- BOLA

Cover Thread 상세 규격

Model	사용 가능 Tubing Size (Max.)	Thread
B318-40	2 x 8.5mm / 1 x 10mm	2 x GL14 / 1 x GL18
B318-64	2 x 10mm / 1 x 16mm	2 x GL18 / 1 x GL25
B318-80	3 x 16mm	3 x GL25



PFA Reaction Vessel PFA 테프론 평저 반응기 하부

Cat. No.	Model	용량	외부 Size	Flange	EA Price ₩
R03-460-051	B320-01	2.4 L	∅ 150 x h 150mm	NW 170 / 146	2,274,000

- 재질 : PFA, -200 °C to 200 °C 내열
- Flat Flange Reaction Vessel
- BOLA



PTFE Reaction Vessel PTFE 테프론 환저 반응기 하부

Cat. No.	Model	용량	외부 Size	Flange	EA Price ₩
R03-460-011	B281-03	500ml	∅ 110 x h 120mm	NW 100	1,130,800
R03-460-015	B281-06	1 L	∅ 110 x h 205mm		1,940,100
R03-460-019	B281-09	2 L	∅ 140 x h 270mm		7,308,100
R03-460-023	B281-12	4 L	∅ 200 x h 290mm	NW 150	11,160,000
R03-460-027	B281-15	6 L	∅ 215 x h 320mm		13,948,400

- 재질 : PTFE, -200 °C to 200 °C 내열
- 테프론 재질로 최상의 내화학성을 가짐
- Flat Flange Reaction Vessel
- BOLA



PTFE 4-N Reaction Vessel Cover PTFE 테프론 4구 반응조 상부

Cat. No.	Model	용량	Joint	EA Price ₩
R04-460-101	B283-02	NW 100	29 / 32	1,731,300
R04-460-105	B283-04	NW 150		2,823,600
R04-460-109	B322-01	NW 170 / 146		1,334,600

- 재질 : PFA, -200 °C to 250 °C 내열
- Side Neck은 15 °의 각도를 가지고 있음
- Flat Flange Lid
- BOLA

