


Flat Flange Gasket 반응기용 가스켓

Cat. No.	Model	Flange	Material	EA Price ₩
R04-460-051	B282-02	NW 100	PTFE	112,800
R04-460-055	B282-04	NW 150		228,700
R04-460-059	B321-01	NW 170 /146	Silicon (Thickness 0.5mm)	130,800

- 재질 : PTFE, Silicon, -60 °C to 230 °C 내열
- 오직 PTFE 재질만 시료에 접촉 되어짐
- BOLA



Flat Flange Joining Clamp 반응기용 클램프

Cat. No.	Model	Flange	Material	EA Price ₩
R04-460-011	B284-02	NW 100	Stainless Steel ①	1,492,900
R04-460-015	B284-04	NW 150		1,293,500
R04-460-019	B323-01	NW 170 /146	Aluminium ②	968,200

- 재질 : Stainless Steel
- 반응기 하부와 상부 결합용 클램프
- BOLA

