


Small Scale Jacketed Reactor, Bench-Top 자켓식 소형 반응 시스템, 300 ml, 500 ml, 1 L

Cat. No.	Model	용량	Velp OHS 100	Manifold Ass.	EA Price ₩
R08-132-001	1949-X-303	300 ml	X	X	18,244,000
R08-132-005	1949-X-300		O		20,040,800
R08-132-009	1949-X-303M		X	O	22,818,500
R08-132-013	1949-X-300M		O		25,814,300
R08-132-017	1949-X-505	500 ml	X	X	17,804,700
R08-132-021	1949-X-500		O		20,198,800
R08-132-025	1949-X-505M		X	O	23,009,200
R08-132-029	1949-X-500M		O		25,917,900
R08-132-033	1949-X-1D	1 L	X	X	17,523,600
R08-132-037	1949-X-1L		O		20,379,200
R08-132-041	1949-X-1DM		X	O	23,242,400
R08-132-045	1949-X-1LM		O		26,151,300

- ChemRxnHub Process System
- 기본구성 : Bench top Reactor kit, Complete with Glassware and Stand
- Manifold (용매제 압력조절 및 입/출구 분리 편리성) 및 모터의 종류에 따라 구분
- OHS 100 Advance Wifi (6b)가 안되는 제품은 Digital Overhead Stirrer (6a)제품으로 구성됨
- Benchtop Hood 안에 설치 가능
- Stand Size : w390 x d440 x h1120 mm
- Chemglass

Standard Component, 구성품

Cat. No.	Model	Description	EA Price ₩
R02-132-544	2091-04	Stirrer Paddle for 10 mm O.D., 50 mm, 60°, PTFE ①	460,200
R02-132-548	2078-05	300 ml & 500 ml용 Stirring Shaft, 10 mm, Polished ②	223,500
R02-132-552	2078-07	1 L용 Stirring Shaft, 10 mm, Polished ②	243,100
R02-132-556	2081-A-04	Anchor Style Stir Paddle, 81 mm x 35 mm, PTFE ③	284,100
R02-132-566	2077-G-21	24/40 Stirrer Bearing, PTFE ④	722,200
R02-132-570	2044-C-14	Tru-Stir Shaft Coupling ⑤	495,300
R02-132-572	2024-10	Digital Overhead Stirrer, 230V (6a)	3,136,300
R02-132-579	2054-100	Velp OHS 100 Advance Stirrer, 230V (6b)	6,037,500
R02-132-137	1929-X11	300 ml Reaction Vessel, Jacketed ⑦	2,320,000
R02-132-141	1929-X13	500 ml Reaction Vessel, Jacketed ⑦	2,424,300
R02-132-145	1929-X15	1000 ml Reaction Vessel, Jacketed ⑦	2,703,600
R02-132-C27	1946-D-01	Lid, Reaction Vessel, 5-neck, 100 mm ⑧	758,800
R02-132-576	141-T-12	100 mm Quick Release Clamp ⑨	423,500
R02-132-580	147-76	O-ring, FEP-encased Stainless	312,800
R02-132-562	1213-J-01	Condenser, Coil, Reflux ⑩	727,300
R02-132-567	1064-01	Adapter, Inert Gas, 24/40, 2-Hose Conn. 180° ⑪	55,700
R02-132-575	1971-M-09	PTFE jacketed S.S. Baffle, 360 mm, 300/500 ml용 ⑫	565,400
R02-132-573	1971-M-10	PTFE Jacketed Baffle 410mm S.S., 1 L용 ⑫	565,400
R02-132-578	1971-M-20	24/40 Micro Baffle Adapter, PTFE ⑬	234,300
R02-132-584	3002-L-24	Stopper, PTFE, 24/40 with Loosening Ring ⑭	288,900
R02-132-526	1949-X-10	Reactor Support Frame and LDPE Drip Tray ⑮	7,793,000
R02-132-586	1947-B100	S.S. Reactor Clamp, 100 mm (16)	659,500
R02-132-598	1969-V-10	30 mm Flat Flange S.S Jacket Adapter, M16 (Only Manifold System)	582,000
R02-132-691	1969-X-10	Inlet + Outlet Manifold Assembly set (17), (18)	8,007,200



Tool-Free Jacket fittings with hose hanger incorporated into CHEMRxnHUB



Unique clamp and support system provides simple installation and quick change-out