

Pilot Scale Process Reactor, Jacketed 자켓식 대형 반응 시스템, 30 L & 50 L

Cat. No.	Model	용량	Motor	Manifold Ass.	EA Price ₩
R02-132-068	1968-630	30 L	Electric	X	58,830,800
R02-132-074	1968-630M			O	64,948,400
R02-132-076	1968-650	50 L		X	61,623,100
R02-132-081	1968-650M			O	68,273,700

- 높이 : Electric 2150 mm (30 L), 2300 mm (50 L) 모터 포함 높이임
- Valve 하단에서 바닥까지의 높이 : 470 mm (30 L), 410 mm (50 L)
- Manifold (용매제 압력조절 및 입/출구 분리 편리성) 따라 구분
- 요청에 의해 Air Type Motor 공급가능
- Chemglass

Standard Component, 구성품

Cat. No.	Model	Description	EA Price ₩
R02-132-792	1215-C-10	High Efficiency Condenser, 45/50, ①	2,385,700
R02-132-796	2097-31	25.4 mm Stirrer Shaft, h1080mm OAL, ②	1,312,700
R02-132-800	2097-53	25.4 mm Stirrer Shaft, h1330mm OAL, ②	1,422,300
R02-132-804	2077-H-06	25.4 mm Teflon Stir Bearing, ③	1,107,600
R02-132-808	149-06	80 mm Glass Cap, ④	326,000
R02-132-812	141-T-10	80 mm S.S. Clamp, Teflon Coated, ⑤	359,300
R02-132-816	147-79	80 mm FEP Encased S.S. Spring O-Ring, ⑥	240,700
R02-132-820	1968-G-01	300 mm Envelope Gasket, ⑦	603,700
R02-132-244	1968-61	30 L Jacketed Reaction Vessel, ⑧	20,650,600
R02-132-252	1968-81	50 L Jacketed Reaction Vessel, ⑧	23,960,200
R02-132-692	1969-A-05	S.S. Adapter, 1" Beaded Pipe to M16x1 Male (2ea), ⑨	199,600
R02-132-692C	HN10HM16PF	1969-A-05 용 Hose Nipple, 10mm Hose-F & 16mmPF	48,200
R02-132-696	1968-66	1" Beaded Pipe Coupling (2ea), ⑩	302,700
R02-132-410	1968-67	2" Beaded Pipe Coupling, ⑪	393,300
R02-132-709	2033-B-25	1/4 HP Brushless DC Motor Only 24v, 3/4" Shaft Diameter ⑫	3,814,100
R02-132-840	2046-90	25.4 mm Shaft Coupling, (Electric), ⑬	814,000
R02-132-500	2033-V-75	Digital Controller Only for Brushless DC Motor, ⑭	4,412,800
R02-132-504	3498-03	1/2" OD Rod, 90° Bend, ⑮	190,700
R02-132-852	9253-20	Universal Kwik Mount Klamp (2ea), ⑯	259,800
R02-132-744	1042-E-10	Thermocouple Adapter, 45/50, ⑰	164,200
R02-132-C72	1968-A-01	7-Neck Lid, 300 mm, ⑱	7,721,600
R02-132-860	2093-21	PTFE Flake Catching Cup, 25.4mm, ⑲	603,700
R02-132-865	1978-P932	1/4" Diameter x 813 mm Long PT100 Temp. Sensor, (20)	1,180,400
R02-132-869	1978-P945	1/4" Diameter x 1140 mm Long PT100 Temp. Sensor, "T", (20)	1,299,800
R02-132-872	2095-23	203 mm OD PTFE Agitator, 60°, (21)	1,236,100
R02-132-880	1968-X-12	New Open Support Frame, 1910 mm OAH, (22)	15,551,000
R02-132-884	1968-X-14	New Open Support Frame, 2120 mm OAH, (22)	15,845,400
R02-132-781	1968-GL-040	Drain Valve Zero Dead Space 2inch BP, GL-40 Thread, (23)	2,115,900
R02-132-888	2048-03	25.4 mm Safety Collar with Teflon Gasket	174,300
R02-132-784	1963-B-10	Beam Torque Wrench	221,900
R02-132-788	1963-B-14	7/16" Deep Socket, 6 Point, 1/4" Drive	15,900
R02-132-901	CG-2097-58	50L-S용 25.4mm Polished Stirrer Shaft 1100mm	1,314,500
R02-132-910	CG-1968-T-50	Large Zero Dead Space Drain Valve	2,571,700
R02-132-833	1969-M-06	Inlet + Outlet Manifold Assembly set	8,201,100

- 시스템 용량에 따라 Reaction Vessel과 Stirring Shaft가 변경되어짐

