

## I Large Process Filter Reactor System, Jacketed 자켓식 대형 반응 시스템, 필터부착형

Cat. No.	Model	용량	Motor	전체높이	EA Price ₩
R02-132-304	1959-F610	10 L	Electric	1830 mm	48,171,300
R02-132-308	1958-F611		Air	1750 mm	43,337,600
R02-132-312	1959-F620	20 L	Electric	1830 mm	49,723,800
R02-132-316	1958-F621		Air	1750 mm	44,920,000

- Complete System Auto-center Stand와 Air 또는 Electronic 모터를 포함
- 200 mm PTFE Filter, 30 um (Pore)
- 높이는 Motor를 포함한 전체 높이임
- Chemglass

### Standard Component, 구성품

Cat. No.	Model	Description	EA Price ₩
R02-132-654	2097-16	19mm Stirrer Shaft, h830mm ①	799,600
R02-132-347	2077-M-03	19 mm PTFE Stir Bearing ②	5,281,400
R02-132-811	2077-M-32	45/50 Lower Joint Adapter ②	1,171,800
R02-132-664	149-01	60 mm Glass Cap ③	278,600
R02-132-668	141-T-11	60 mm S.S. Clamp, PTFE Coated, ④	319,300
R02-132-672	147-80	60 mm Perfluoro O-Ring, ⑤	656,900
R02-132-676	147-84	200 mm Perfluoro O-Ring (2ea) ⑥	1,490,700
R02-132-682	1959-10	10L Jacketed Filter Reaction Vessel ⑦	7,927,300
R02-132-690	1959-20	20L Jacketed Filter Reaction Vessel ⑦	9,720,600
R02-132-692	1969-A-05	S.S. Adapter, 1" Beaded Pipe to M16×1Male (2ea), ⑧	199,600
R02-132-692C	HN10HM16PF	1969-A-05 용 Hose Nipple, 10mm Hose-F & 16mm PF	48,200
R02-132-696	1968-66	1" Beaded Pipe Coupling (2ea) ⑨	302,700
R02-132-709	2033-B-25	1/4 HP Brushless DC Motor Only 24v, 3/4" Shaft Diameter, ⑩	3,814,100
R02-132-708	2025-20	Air Motor, 3/8HP, ⑩	2,277,800
R02-132-712	2025-10	Filter-Regulator-Lubricator (Air Motor Only)	695,500
R02-132-518	2044-30	19 mm Tru-Stir® Shaft Coupling, (Electric), ⑪	575,500
R02-132-720	2044-16	19 mm Tru-Stir® Shaft Coupling, (Air), ⑪	575,500
R02-132-500	2033-V-75	Digital Controller Only for Brushless DC Motor, 120v (Electric Only) ⑫	4,412,800
R02-132-504	3498-03	1/2" OD Rod, 90° Bend ⑬	190,700
R02-132-512	9253-10	Universal Kwik Mount Klamp ⑭	148,400
R02-132-744	1042-E-10	Thermocouple Adapter, 45/50 ⑮	164,200
R02-132-C68	1946-50	5-Neck Lid, 200 mm (16)	3,014,300
R02-132-752	141-T-15	200 mm S.S. Clamp, PTFE Coated, (17)	573,500
R02-132-763	11978-P922	10L용 1/4" Diameter×560 mm Long PT100 Temp. Sensor (18)	1,106,500
R02-132-761	1978-P926	20L용 1/4" Diameter×660 mm Long PT100 Temp. Sensor (18)	1,107,200
R02-132-764	2095-05	127 mm OD PTFE Agitator, 60° (19)	1,187,700
R02-132-766	2096-F-08	127 mm OD PTFE Agitator, 90° (20)	925,900
R02-132-754	1959-F200	Filter Support Assembly, 200 mm (21)	8,264,700
R02-132-996	1959-F203	PTFE Filter 30um, 200 mm	961,300
R02-132-776	1965-X-150	Auto-Center Support Frame, 1610 mm OAH (22)	10,267,500
R08-132-836	1964-A-10	1" Beaded Pipe to 1/2" NPT M, Glass Filled PTFE (23)	205,900
R02-132-782	1959-A-01	Drain Valve, 0-14, 1/2" Side Arm (24)	528,600
R02-132-784	1963-B-10	Beam Torque Wrench	221,900
R02-132-788	1963-B-14	7/16" Deep Socket, 6 Point, 1/4" Drive	15,900

- 시스템 용량에 따라 Reaction Vessel과 Stirring Shaft가 변경되어짐

