

Small Scale Filter Reactor System, Unjacketed 소형 반응 시스템, 필터부착형

Cat. No.	Model	용량	Motor	EA Price ₩
R08-132-200	1959-U1000	1 L	Electric	30,436,400
R08-132-204	1959-U5000	5 L		34,645,700

- Complete System : Stand와 Electronic 모터 포함
- 100 mm PTFE Filter, 30 um (Pore)
- Benchtop Hood 안에 설치 가능
- Chemglass

Standard Component, 1 L 반응기 구성품

Cat. No.	Model	Description	EA Price ₩
R02-132-800	2095-X-102	50 mm PTFE Agitator, 60° Impeller Angle ①	1,422,300
R08-132-804	2075-U10	10 mm Stirring Shaft, 480 mm OAL, 4 Blade ②	485,200
R02-132-566	2077-G-21	10 mm PTFE Stirrer Bearing 24/40 ③	722,200
R02-132-570	2044-C-14	Compact Stirrer Shaft Coupling ④	495,300
R02-132-579	2054-100	Overhead Stirrer Remote 220V ⑤	6,037,500
R08-132-816	1959-U02	1 L UnJacketed Filtration Vessel ⑥	736,700
R08-132-820	1946-B-01	Lid, 5-Neck, 100 mm Flange, 24/40 ⑦	785,100
R02-132-576	141-T-12	Quick Release Clamp, 100 mm Flange, PTFE Coated ⑧	423,500
R08-132-824	147-82	100 mm Perfluoro O-ring	771,600
R08-132-828	1949-F100	100 mm Filter Support Assembly, w/o Filter ⑨	3,757,000
R08-132-832	1949-F103	100 mm PTFE Filter, 30 um ⑩	301,400
R08-132-836	1964-A-10	Glass-Filled PTFE Drain Adapter ⑪	205,900
R02-132-782	1959-A-01	Drain Valve ⑫	528,600
R02-132-696	1968-66	1" Beaded Pipe Coupling ⑬	302,700
R02-132-584	3002-L-24	PTFE Stopper, 24/40 w/Loosening Ring ⑭	288,900
R08-132-844	3059-F100	100 mm Filter Support Assembly Jack ⑮	1,951,300
R02-132-526	1949-X-10	Reactor Support Frame (16)	7,793,000
R02-132-586	1947-B100	100 mm SUS Reactor Clamp (17)	659,500

Standard Component, 5 L 반응기 구성품

Cat. No.	Model	Description	EA Price ₩
R02-132-632	2095-04	4" PTFE Agitator, 60° Impeller Angle ①	1,152,500
R08-132-806	2097-U50	19 mm Stirrer Shaft, 4 Blade Glass ②	924,500
R02-132-616	2077-G-24	19 mm PTFE Stirrer Bearing, 45/50 ③	1,394,600
R02-132-571	2044-C-16	Compact Stirrer Shaft Coupling ④	523,400
R02-132-579	2054-100	Overhead Stirrer Remote 220V ⑤	6,037,500
R08-132-818	1959-U05	5 L UnJacketed Filtration Vessel ⑥	1,203,700
R08-132-825	1946-B-11	5-Neck Lid, 150mm Flange, 45/50 ⑦	1,073,500
R02-132-624	141-T-14	150 mm S.S Clamp, PTFE Coated ⑧	505,900
R08-132-826	147-83	150 mm Perfluoro O-ring	1,309,600
R08-132-830	1949-F150	150 mm Filter Support Assembly, w/o Filter ⑨	4,883,700
R08-132-834	1949-F153	150 mm PTFE Filter, 30um ⑩	539,800
R08-132-836	1964-A-10	Glass-Filled PTFE Drain Adapter ⑪	205,900
R02-132-782	1959-A-01	Drain Valve ⑫	528,600
R02-132-696	1968-66	1" Beaded Pipe Coupling ⑬	302,700
R02-132-584	3002-L-24	PTFE Stopper, 24/40 w/Loosening Ring ⑭	288,900
R08-132-846	3059-F150	150 mm Filter Support Assembly Jack ⑮	2,019,000
R02-132-526	1949-X-10	Reactor Support Frame (16)	7,793,000
R02-132-587	1947-B150	150 mm SUS Reactor Clamp (17)	696,100

