

I Large Process Filter Reactor System, Unjacketed 대형 반응 시스템, 필터부착형

Cat. No.	Model	용량	Motor	전체높이	EA Price ₩
R02-132-334	1959-F810	10 L	Electric	1830 mm	42,701,700
R02-132-338	1958-F811		Air	1750 mm	38,469,000
R02-132-342	1959-F820	20 L	Electric	1830 mm	44,885,700
R02-132-346	1958-F821		Air	1750 mm	40,551,900

- Complete System Auto-center Stand와 Air 또는 Electronic 모터를 포함
- 200 mm PTFE Filter, 30 um (Pore)
- 높이는 Motor를 포함한 전체 높이임
- Chemglass

Standard Component, 구성품

Cat. No.	Model	Description	EA Price ₩
R02-132-654	2097-16	19mm Stirrer Shaft, h830mm ①	799,600
R02-132-347	2077-M-03	19 mm PTFE Stir Bearing ②	5,281,400
R02-132-811	2077-M-32	45/50 Lower Joint Adapter ②	1,171,800
R02-132-664	149-01	60 mm Glass Cap ③	278,600
R02-132-668	141-T-11	60 mm S.S. Clamp, PTFE Coated ④	319,300
R02-132-672	147-80	60 mm Perfluoro O-Ring, ⑤	656,900
R02-132-676	147-84	200 mm Perfluoro O-Ring (2ea) ⑥	1,490,700
R02-132-785	1959-U10	10L Unjacketed Filter Reaction Vessel ⑦	2,548,300
R02-132-790	1959-U20	20L Unjacketed Filter Reaction Vessel ⑦	5,001,900
R02-132-696	1968-66	1" Beaded Pipe Coupling (2ea) ⑧	302,700
R02-132-709	2033-B-25	1/4 HP Brushless DC Motor Only 24v, 3/4" Shaft Diameter	3,814,100
R02-132-708	2025-20	Air Motor, 3/8HP ⑨	2,277,800
R02-132-712	2025-10	Filter-Regulator-Lubricator (Air Motor Only) ⑨	695,500
R02-132-518	2044-30	19 mm Tru-Stir® Shaft Coupling, (Electric) ⑩	575,500
R02-132-720	2044-16	19 mm Tru-Stir® Shaft Coupling, (Air) ⑩	575,500
R02-132-508	3499-02	3 m Coiled Extension Cord, "T" ⑪	99,800
R02-132-500	2033-V-75	Digital Controller Only for Brushless DC Motor (Electric Only) ⑫	4,412,800
R02-132-504	3498-03	1/2" OD Rod, 90° Bend ⑬	190,700
R02-132-512	9253-10	Universal Kwik Mount Klamp ⑭	148,400
R02-132-744	1042-E-10	Thermocouple Adapter, 45/50 ⑮	164,200
R02-132-C68	1946-50	5-Neck Lid, 200 mm (16)	3,014,300
R02-132-752	141-T-15	200 mm S.S. Clamp, PTFE Coated, (17)	573,500
R02-132-763	1978-P922	1/4" Diameter×560 mm Long PT100 Temp. Sensor (18)	1,106,500
R02-132-761	1978-P926	1/4" Diameter×660 mm Long PT100 Temp. Sensor (18)	1,107,200
R02-132-764	2095-05	127 mm OD PTFE Agitator, 60° (19)	1,187,700
R02-132-766	2096-F-08	127 mm OD PTFE Agitator, 90° (20)	925,900
R02-132-754	1959-F200	Filter Support Assembly, 200 mm (21)	8,264,700
R02-132-996	1959-F203	PTFE Filter 30um, 200 mm	961,300
R02-132-776	1965-X-150	Auto-Center Support Frame, 1610 mm OAH (22)	10,267,500
R08-132-836	1964-A-10	1" Beaded Pipe to 1/2" NPT M, Glass Filled PTFE (23)	205,900
R02-132-782	1959-A-01	Drain Valve, 0-14, 1/2" Side Arm (24)	528,600
R02-132-784	1963-B-10	Beam Torque Wrench	221,900
R02-132-788	1963-B-14	7/16" Deep Socket, 6 Point, 1/4" Drive	15,900

- 시스템 용량에 따라 Reaction Vessel과 Stirring Shaft가 변경되어짐

